
CAT

Ingenieurbüro CAT
M. Zipperer GmbH

Instruction Manual

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Content

1. User Instructions	6
1.1. Important Instructions for your safety.....	6
1.2. Danger symbols in this operating manual.....	6
Danger symbols:	6
Danger levels	7
1.3. Warning signs on the device	7
2. General safety warnings and instructions	8
3. Scope of Delivery	10
4. General Information	10
5. Intended Use.....	10
6. Setting up and Starting up the Instrument.....	11
6.1. Unpacking the Instrument.....	11
6.2. Setting up the Instrument	11
6.3. Electrical Connection.....	12
6.4. Rear Connections	12
6.5. Connection of the Pt100 Probe.....	13
6.6. Probe Stand	13
6.7. Extension Connectors.....	14
7. The Control Panel.....	14
7.1. Description of the Control Panel Functions	14
7.2. Switching the Device On and Off.....	16
8. Heating.....	16
8.1. Temperature Control with external Pt100 Probe	16
8.2. Temperature Control without Pt100 Probe (Hot Plate Control).....	17
8.3. Plate Limit (Heating with Pt100 probe).....	18
8.4. Safety Temperature.....	18
9. Fahrenheit-Celsius Setting	19
10. Stirring.....	20
11. Timer Controlled Operation.....	20
12. Additional Safety Functions	20
12.1. Differential Alarm	21
12.2. Out of Liquid Check	21
12.3. Failure of External Pt100 Probe	21
12.4. Internal Over Heating Monitoring.....	21
13. Error Messages.....	22
13.1. Other Messages.....	23

14.Serial RS485 Interface	23
14.1. Pin assignment of the RS485 Interface	23
14.2. Connection between PC and the device unit.....	24
14.3. Settings of the RS485 Interface	24
14.4. Addressing of the device unit.....	25
Addressing via front panel.....	25
Addressing via RS485 command.....	25
14.5. RS485 Data transfer.....	26
14.6. Format of a RS485 command.....	26
14.7. Format of a RS485 Handshake	26
14.8. RS485 Commands.....	27
14.9. RS485 Examples of use.....	29
Example 1: Setting of values.....	29
Example 2: Reading of actual values	30
15.Maintenance and Cleaning	31
16.Warranty and Liability.....	31
17.Dismantling, Transport and Storage	32
17.1. Dismantling.....	32
17.2. Transport and Storage.....	32
18.Technical Data.....	32
19.Repair Return Form.....	35
37.Declaration of Conformity / Konformitätserklärung.....	64

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1. User Instructions

1.1. Important Instructions for your safety

- Every user must read and understand this manual completely before use. Failure to do so can result in serious injury or death.
- Comply with all safety and accident-prevention regulations applicable to laboratory work.
- Follow general instructions for hazard prevention and general safety instructions, e.g. wear protection clothing, eye protection and gloves.
- This operating manual is part of the product. Thus, it must always be easily accessible.
- This instruction sheet does not purport to address all of the safety problems which might result from the use of this device, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult their authorized safety advisors and establish appropriate health and safety practices and then determine the application of regulatory limitations prior to use.
- Enclose this operating manual when transferring the device to another place.
- If this manual is lost, please request another one. Please contact your dealer or

Ingenieurbüro CAT
M. Zipperer GmbH
Wettelbrunner Str. 6
D-79282 Ballrechten-Dottingen
Tel.: ++49-(0)7634-5056-800
Fax: ++49-(0)7634-5056-801
www.cat-ing.de
info@cat-ing.de

1.2. Danger symbols in this operating manual

The safety instructions in this manual appear with the following danger symbols and danger levels:

Danger symbols:

	Hazard point		Electrical shock
	Risk of fire		Explosion

	Crushing		Hot surface
	Bio hazard		Magnetic field
	Chemical hazard		Material damage

Danger levels

	Will lead to severe injuries or death
	May lead to severe injuries or death
	May lead to light to moderate injuries
	May lead to material damage

1.3. Warning signs on the device

	<p> WARNING</p> <p>This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.</p>
	<p> WARNING</p> <p>Beware of the effect of the magnetic field on cardiac pacemakers and data storage media.</p>
	<p> WARNING</p> <p>This symbol indicates that there are hot surfaces on the device. It refers to the hotplate of the device which can reach up to 500°C during operation causing a hazard resulting in serious burns.</p> <ul style="list-style-type: none"> • Avoid contact with the hotplate, vessel and heating medium. • Avoid storage of flammable substances near the instrument.

2. General safety warnings and instructions

	<p>⚠ DANGER Risk of explosion.</p> <ul style="list-style-type: none"> • Do not operate the device in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof. • Do not use this device for processing any substances which could generate an explosive atmosphere. • Do not use this device to process any explosive or highly reactive substances. • Do not use flammable or explosive substances near the instrument.
	<p>⚠ DANGER Electric shock as a result of penetration of liquid.</p> <ul style="list-style-type: none"> • Do not allow any liquids to penetrate the inside of the housing • Switch off the device and disconnect the power plug before starting cleaning or disinfection work. The On/Off Switch on the device does not disconnect the device from the power source. • Use only closed tubes. • Only plug the device back in if it is completely dry, both inside and outside.
	<p>⚠ WARNING Risk from incorrect supply voltage</p> <ul style="list-style-type: none"> • Only connect the device to a AC power source with an protective earth (PE). • Only connect the device to voltage sources which correspondent to the electrical requirements on the type label.
	<p>⚠ WARNING Electric shock due to damage to device or mains cable</p> <ul style="list-style-type: none"> • Only connect the device to the mains supply if the device and the mains cable are undamaged • Only use devices that have been properly installed or repaired. • In case of danger, disconnect the device from the mains supply by pulling the power plug from the mains socket or by using the isolating device intended for this purpose (e.g. emergency stop switch)
	<p>⚠ WARNING Lethal voltage inside the device</p> <ul style="list-style-type: none"> • Do not open the device. • Ensure that the housing is always closed and undamaged so that no parts inside the housing can be contacted by accident. • The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely. • Do not allow any liquids to penetrate the inside of the housing. • Repairs are only to be carried out by trained service technicians.
	<p>⚠ WARNING Damage to health due to corrosive or aggressive chemicals</p> <ul style="list-style-type: none"> • Observe all markings on the reagent bottles. • Always check the instrument for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation. <p>Leaking solutions may endanger persons and materials</p>

	<ul style="list-style-type: none"> ● Use proper connecting vessels, protective clothing and gloves. ● Avoid splashes ● When dispensing, maintain a physical distance between the instrument and the body. ● Dangerous and fuming chemicals must be dispensed in a fume hood. ● Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction. ● Always use the instrument in such a manner that neither the operator, nor any other person is endangered.
	<p>▲WARNING Damages to health due to infectious liquids and pathogenic germs.</p> <ul style="list-style-type: none"> ● When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets and the manufacturer's application notes. ● Wear personal protective equipment ● For comprehensive regulations about handling germs or biological material of the risk group II or higher, please refer to the "Laboratory Biosafety Manual" in its respectively current valid version from the World Health Organisation
	<p>▲WARNING Damage to health due to contaminated device and accessories</p> <p>In the following cases, sample material can be released:</p> <ul style="list-style-type: none"> - improperly sealed tubes - unstable tubes - high vapour pressure of the content so that the seal of the tubes can spring open - damaged sealings - smashed glass tubes <ul style="list-style-type: none"> ● Only mix in closed tubes ● Observe the nationally prescribed safety environment when working with hazardous, toxic and pathogenic samples. Pay particular attention to personal protective equipment (gloves, clothing, goggles, etc.), extraction, and the safety class of the lab. ● Decontaminate the device and the accessories before storage and shipping.
	<p>▲WARNING Risk of fire</p> <ul style="list-style-type: none"> ● Do not use this device to process any highly flammable liquids
	<p>▲CAUTION Poor safety due to inadequate fixing of the unit</p> <ul style="list-style-type: none"> ● Ensure that the unit is firmly attached to a solid stand.
	<p>▲CAUTION Poor safety due to incorrect accessories and spare parts.</p> <p>The use of accessories and spare parts other than recommended by Ingenieurbüro CAT, M. Zipperer GmbH may impair the safety, function and precision of the device. Ingenieurbüro CAT, M. Zipperer GmbH cannot be held liable or accept any liability</p>

	<p>for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.</p> <ul style="list-style-type: none"> ● Only use accessories and spare parts recommended by Ingenieurbüro CAT, M. Zipperer GmbH
	<p>CAUTION Crush hazard due to moving parts</p> <ul style="list-style-type: none"> ● Do not replace any consumables as long as the device is running. ● Do not open the coverage as long as the device is running

3. Scope of Delivery

Unpack the delivered contents carefully and check to see that nothing appears physically damaged or is missing. Your delivery should contain the following:

- | | | |
|-------------------------------------|--------------------------|--------------------------|
| - 1 Magnetic Hotplate Stirrer M 21, | PN: 60361-0000 (230V) or | PN: 60361-0001 (115V) or |
| - 1 Magnetic Hotplate Stirrer M 22, | PN: 60362-0000 (230V) or | PN: 60362-0001 (115V) or |
| - 1 Magnetic Hotplate Stirrer M 23, | PN: 60363-0000 (230V) or | PN: 60363-0001 (115V) |
| - 1 Instruction Manual | | |

Optional:

- | | |
|--|----------------|
| - 1 Pt100 Temperature Probe, | PN: 60278-0000 |
| - Stand, | PN: 60663-0000 |
| - Probe Clamp, | PN: 8B00562400 |
| - Clamp, | PN: 60668-0000 |
| - RS485-RS232 Adapter for RS485 Communication, | PN: 61703-0000 |
| - RS 232 to USB Adapter, | PN: 30244-0001 |
| - 9 pol. Sub-D extension cord, | PN: 30275-0051 |

4. General Information

Our Magnetic Hotplate Stirrers are designed according to protection class 1. They were manufactured and tested according to DIN EN 61010.

According to these regulations, the units are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

5. Intended Use

These Magnetic Hotplate Stirrers are mainly to be used to heat and stir hydrous solutions and oil not containing inflammable parts in glass vessels. They are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

6. Setting up and Starting up the Instrument

6.1. Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

If the instrument is not damaged and all parts are complete you may start to operate the device after reading the instruction manual.

6.2. Setting up the Instrument

Please place the device on a fire-proof horizontal and even surface. A minimum distance of 50cm to inflammable materials should always be observed.

	▲ DANGER The device may not be operated in explosion-prone areas.
	▲ WARNING The power cable may not come in contact with the hotplate.
	▲ WARNING The device is not to be used without supervision.
	▲ WARNING When connecting the instrument to an AC power outlet ensure that your local supply voltage matches the indication on the instrument.
	▲ WARNING Beware of the effect of the magnetic field on cardiac pacemakers and data storage media.

6.3. Electrical Connection


- A earthing-pin plug (DIN 49441 CEE 7/VII10/ 16 A 250 V, a standard plug in Germany, Austria, the Netherlands, Belgium, France, Norway, Sweden, Finland, Denmark, Portugal and Spain) is standard on all instruments. For North America instruments feature standard US plugs (NEMA Pub.No.WDI1961 ASA C 73.1. 1961 page 8 15A 125V),, for the UK with a standard UK plug BS 1363.
- When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.
- The instrument is earthed as supplied. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!
- When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

	⚠WARNING Before connection the instrument to the mains ensure that the rocker switch (1) is set to the „OFF“- position.
--	---

Connect the mains Power Cable to the mains power connector. Ensure that your local supply voltage matches the indication on the instrument

6.4. Rear Connections


(Rear view of the device)

6.5. Connection of the Pt100 Probe

A DIN Connector on the rear of the instrument is intended for the connection of a Pt100 temperature Probe (Order No. 60278-0000). Pt1000 Probes are not suitable.

Verbindung eines PT100 Sensors:


(Connection with a Pt100 probe)

	CAUTION Ensure that the cable of the Pt100 probe does not come in contact with the hotplate.
	NOTICE Use Pt100 probes with electrically isolated tips only. Probes that allow for an electrical connection to the metal chassis of the device may distort temperature measurement.
	NOTICE Ensure that the Pt100 probe tip is submersed to a depth of at least 50mm in the heating medium to allow accurate measurement.

6.6. Probe Stand

Use the M6 Thread on the unit for the optional probe stand (Order No. 60663-0000).

(Mounted probe stand)


6.7. Extension Connectors

These are intended for the connection of external modules or RS485-RS232 adapters for RS485 communication to a Personal Computer. Please contact your supplier or manufacturer for more information and order numbers.

	NOTICE Never connect the RJ45 jack with a network jack of a PC, router, switch or something like that. Only connect original accessories to the "RS485 In" jack! Any other use may result in damage of the device or damage to third party devices!
---	--

7. The Control Panel

The control wheel on the left is used for changing set values. Any values and messages to be displayed are shown on the 4 digit 7-segment display. The turning knob on the right is used for setting the speed of the stirring motor. When switching on the unit the set motor speed will be displayed on the LED.

(Control panel of the device)

7.1. Description of the Control Panel Functions

7-Segment LED Display

This LED display is used for showing different values and messages (temperature, motor speed, timer and error messages)

Button: „Plate On/Off“

To switch on the heating plate after entering the desired set-values this button needs to be pressed. Switching the plate off is also done by pressing this button.

Button: „View RPM“
Show the actual set motor speed

Button: „Plate Temp“
Display the actual plate temperature or press again for setting the set point temperature. When used with Pt100 probe: display the actual plate temperature or press again for setting the plate limit temperature.


Button: „Safety Temp“
Display or set the safety temperature


Button: „Probe Temp“
Display or set the probe temperature (Pt-100)

Button: „Timer“
Display or set the remaining timer minutes

Motor Speed Wheel
This wheel is used to set the speed of the stirring motor. The set speed is displayed on the LED display while this wheel is turned (opto-electronic speed control)

Encoder Wheel with Push Button
This encoder wheel is used to change different set values. Turning clockwise will increment values, while turning anti-clockwise will decrement. Pressing the encoder wheel push button will confirm a set value. Any error messages displayed at start can also be reset by pressing the encoder wheel push button.

Status LED „Plate On/Off“

This LED is lit when the hotplate is on. It indicates that the unit is trying to reach or keep the entered set temperature for hotplate or probe.

Status LED „Heater“

This LED indicates that the hotplate is currently being exposed to an energy pulse.

Status LED „Hot“

When the hotplate is hotter than 60°C this LED lights up.

Status LED „Timer“

When the timer is activated this LED lights up.

7.2. Switching the Device On and Off

	▲WARNING The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.
---	---

Make sure that you have read the instruction manual carefully before operating the device. It is also necessary to check whether the device is in accordance with the setup instructions above.

	If the device is in standby mode, the device can be switched on by pressing the I/O key. In switched on mode you can press the I/O key again to switch the device back in standby mode.
---	---

8. Heating

The Magnetic Hotplate Stirrer allows 2 different modes of operation:

- **Temperature control with external Pt100 probe (recommended mode)**
- **Temperature control of the hotplate (no Pt100 probe connected)**

For exact temperature control of a liquid, the use of an external Pt100 probe is recommended to regulate the liquid temperature as accurately as possible.

	▲WARNING Always switch off the device after use!
	▲CAUTION As long as the Heater LED is flashing or lights, the hotplate is supplied with energy. The Hot LED lights when the hotplate is hotter than 60°C. You may get burned!
	▲CAUTION After switching off the device the hotplate stays hot! Therefore, after switching off the device, the display will show the message „HOT“ as long as the hot plate temperature exceeds 60°C and the power cord is still plugged in.

8.1. Temperature Control with external Pt100 Probe

	▲WARNING Please ensure that the cable of the Pt100 probe does not come in contact with the hotplate.
---	---

We recommend the use of an external Pt100 probe (Order No. 60278-00) to allow for accurate temperature control of the heating medium.

- **Connect** the Pt100 probe with the DIN connector on the rear of the Magnetic Hotplate Stirrer unit.
- **Place** the Pt100 probe in the vessel (min. 5cm submersion)

Press Power On/Stby to switch on the device. Any error messages displayed at start can be reset by pressing the encoder wheel button.

	When „ A “ (Amount) is shown as the first digit in the display turn the encoder wheel until the desired volume (in litres) of the heating medium is shown. Press the encoder wheel button to confirm the volume in the vessel.
	Press Probe Temp (press again if the P does not blink). Set the desired setpoint temperature value by turning the encoder wheel. Press the <i>Probe Temp</i> key again or press the encoder wheel button to confirm the set value or wait 3 seconds.
	Press Plate On/Off to switch on the hotplate (Status LED above the button will now be lit).
	CAUTION When a Pt100 probe is connected, the value set with „Plate Temp“ does not refer to setpoint of the hotplate but to the plate limit temperature (Plate Limit – See section 8.3). Normally when using a Pt100 probe the plate limit should be set to its maximum value (M 21: 350°C, M 22: 380°C and M 23: 500°C) to ensure optimal temperature control.
	CAUTION When the Pt100 probe is connected to the device, the setpoint temperature of the probe will be set automatically to 0°C, the hotplate will be switched off and the plate limit (see section 8.3) temperature will be set to its maximum value.
	CAUTION Ensure that the cable of the Pt100 probe does not come in contact with the hotplate.

8.2. Temperature Control without Pt100 Probe (Hot Plate Control)

For certain applications it could be necessary to control the hotplate temperature instead the probe temperature. This is only possible if no Pt100 probe is connected to the device.

Press Power On/Standby to switch on the device. Any error messages displayed at start can be reset by pressing the encoder wheel button.

Press Plate Temp (press again if the **H** does not blink). Set the desired setpoint temperature value of the hotplate by turning the encoder wheel. Press the *Plate Temp* key again or press the encoder wheel button or wait 3 seconds to confirm the set value.

Press Plate On/Off to switch on the hotplate (Status LED above the button will now be lit).

CAUTION If a Pt100 sensor is connected, the set value “Plate Temp” doesn’t refer to the set value of the hotplate, but to the temperature limit of the hotplate (see chapter 8.3)
If the Pt100 sensor is being unplugged, then the set temperature of the hotplate is automatically set to 0°C.

8.3. Plate Limit (Heating with Pt100 probe)

	<p>▲WARNING The use of the hotplate temperature limitation results in a longer heat-up time of the heating medium. Therefore in normal operation, the hotplate temperature limit should be set to the maximum value to achieve optimal heat-up times. Setting the hotplate temperature limit to the maximum value allows the microprocessor to freely select the optimum hotplate temperature to reach the set probe temperature as fast as possible. However, when operating with flammable liquids it may be required to limit the maximum allowed hotplate temperature to a safe value.</p>
---	---

For normal operation, the plate limit temperature should be set to the maximum value to achieve optimal heat-up times. Setting the limit temperature to the maximum value allows the microprocessor to freely select the optimum hotplate temperature to reach the set probe-temperature. However, when operating with flammable liquids it may be required to limit the maximum allowed hotplate temperature to a safe value.

Note: Is the hotplate temperature limited to a too low value the set temperature for the fluid cannot be reached or an extremely long heat-up time will be the case.

If the setpoint of the external Pt100 probe is changed, the plate limit temperature will automatically be set to its maximum value. The new value appears in the display and could be changed to the desired value directly by using the encoder wheel.

To set the hotplate limit temperature when Pt100 probe is connected:

Press *Plate Temp* (press again if the **H** does not blink). Set the desired plate limit temperature by turning the encoder wheel. Press the *Plate Temp* key again or press the encoder wheel button or wait 3 seconds to confirm the set value.

8.4. Safety Temperature

	<p>▲WARNING If there is no security risk, the user should set the Safety Temperature to a value as high as possible to avoid an unintentional switch off of the device in normal operation. In safety critical applications the proper use of the safety temperature is advised.</p>
	<p>▲WARNING If a Pt100 probe is connected to the device, the Safety Temp refers to the temperature of the heating medium, otherwise to the hotplate temperature.</p>

A safety temperature can be set to shut down the device if there is an error or to prevent an unexpected heat up of the medium.

The safety temperature is set automatically by the unit when the Pt100 probe or hotplate setpoint temperature is changed. If the unit is operated with an external Pt100 probe and the setpoint is changed the device automatically set the safety temperature 15°C above the new setpoint. If the liquid temperature exceeds the safety temperature (temperature measured with the external probe) the unit will shut down.

Conversely, if the unit is operated without an external probe and the hotplate setpoint is changed the device automatically set the safety temperature 15°C above the new hotplate setpoint. If the temperature measured at the hotplate (internal sensor) exceeds the safety temperature, the unit will also shut down.

If the device automatically changes the safety temperature the new value is shown in the display for a short time. After the device automatically set the new value another safety temperature can be set manually by user.

To set the safety temperature:

Press **Safety Temp** (the  starts to blink). Set the desired Safety temperature by turning the encoder wheel. Press Safety Temp again to confirm the set value or press the encoder wheel button or wait 3 seconds. The set value will be effective until setpoint of *Probe Temp* or setpoint of *Plate Temp* are changed and the *Safety Temp* is automatically set.

	<p>▲WARNING When a Pt100 probe is used the safety temperature refers to the temperature measured by the Pt100 probe. Otherwise the safety</p>
	<p>▲WARNING In safety critical applications the proper use of the safety temperature is advised. Whenever the setpoint temperature of the external probe (Pt100) or the hotplate is changed the current safety temperature is set to 15°C above the setpoint of the Probe Temp or the Plate Temp value. To set the safety temp as close to the actual setpoint value as possible, set the probe or plate temperature first, before changing the safety temperature value.</p>

9. Fahrenheit-Celsius Setting

This Magnetic Hotplate Stirrer can show temperature values in degrees Celsius as well as Fahrenheit. When the device is switched on in Celsius mode the normal start sequence will be displayed. In Fahrenheit Mode the message **FAH** will be displayed after switching on the device.

You can change between Celsius and Fahrenheit mode in the following way:

Press **Power On/Stby** to switch on the device.

Shortly **after** switching the device on **press the Probe Temp Key and hold.**


The display now shows a flashing **U** and you can change between the settings **CEL** (Celsius) or **FAH** (Fahrenheit) by turning the encoder wheel. Press the encoder wheel to confirm the selected temperature mode. The unit will remember the selected temperature mode until it is changed again by the user.

10. Stirring

Make sure that the vessel used is permeable to magnetic field lines (e.g. glass, ceramics or stainless steel)

To set the motor speed:

Set the motor speed by turning the motor speed wheel

The set motor speed is shown on the LED display while the motor speed wheel is being turned. Turn the wheel entirely to the left to switch the motor off.

Note:

Optimal magnetic coupling is acquired with stirring bars with lengths of 30 – 70 mm and diameters of 5 – 10 mm. The appropriate stirring bar for an application must be determined by the user experimentally.

CAUTION

Use a suitable vessel (e.g. Erlenmeyer flask) and beware of high stirring speeds and high liquid levels to avoid splashing. Make sure that suitable protective clothing and eyewear is used.

11. Timer Controlled Operation

This device has an internal timer function. The timer allows for automatic switching off. The switch off time can be selected in the following way.

E.g. a set time of 10 minutes will cause the shut-off of the unit after 10 minutes.

Press **Timer** key (press again if the  does not blink). Turn the **encoder wheel** to select the desired switch off time (remaining on-time in minutes). Press the Timer key again or press the encoder wheel button or wait 3 seconds to confirm.

Minimum time until Switch off: 1 minute

Maximum time until Switch off: 999 minutes = 16h 39min

12. Additional Safety Functions

The following safety functions provide additional protection from uncontrolled heating. If a safety problem is detected, the device turns off automatically. The reason for the shutdown is displayed on the LED-display. When switching on the device next time, you have to confirm the displayed safety problem by pressing the encoder wheel button.

12.1. Differential Alarm

Breakage of the vessel could create a hazardous situation. To monitor this situation problem this instrument has a '*differential alarm*' feature. The '*differential alarm*' switches the unit off if the probe has fallen out of the liquid (e.g. breakage of glass). The '*differential alarm*' operates if the temperature fall detected by the Pt100 probe exceeds the internal differential alarm sensitivity. The unit then switches off.

12.2. Out of Liquid Check

The unit monitors, whether the Pt100 probe tip is immersed in the liquid in the vessel on the hotplate. The microprocessor checks whether the probe temperature changes in relation to the temperature changes of the hotplate. If the microprocessor detects no increase in the liquid temperature over a certain time, although the hotplate temperature is rising, a warning message is shown on the display and the unit shuts down.

Note:

This safety feature is only active if the setpoint temperature of the probe is at least 10°C higher than the actual measured temperature.


You can change the actual Out of Liquid sensitivity of the unit easily via the front panel. For changing the sensitivity you have to do the following steps:

Press **Power On/Standby** to switch on the device.

Shortly **after** switching the device on **press the Plate Temp Key and hold**.

The display now shows a flashing  and the actual Out of Liquid sensitivity. Now you can change the value (0 ... 100, Standard M 21 and M 22: 20, Standard M 23: 15) by turning the encoder wheel. Smaller values decrease the sensitivity. If you select the value 0 the Out of Liquid function is disabled. Press the encoder wheel to confirm the selected sensitivity. The unit will remember the selected value until it is changed again by the user.

12.3. Failure of External Pt100 Probe

Disconnection or failure of the external Pt100 probe after the device is switched on will be detected and cause the unit to switch off.

12.4. Internal Over Heating Monitoring

The temperature inside the device is constantly monitored. A safety mechanism will switch the device off in case the internal temperature exceeds a safe value.

13. Error Messages

This unit has an on-board self-diagnostic program that detects possible failures. In case of an error message the program switches off the unit. To enable the user to prevent the problem from re-occurring, the unit will display any of the following error messages as it is switched on again. To reset any error message press the encoder wheel button:

LED Display:	Error name:	Description/Reasons :
HErrr	hotplate failure	The hotplate temperature exceeded its maximum value The hotplate temperature sensor may be broken
OFFEr	Off Error	The device was not switched off with the power button. Disconnecting the power cable or a mains power failure may be the source for error.
H4AF	Plate > Safety	The measured hotplate temperature was higher than the set safety temperature for the hotplate (Safety Temp.)
P4AF	Probe > Safety	The measured probe plate temperature was higher than the set safety temperature for the Pt100 probe (Safety Temp.)
CLOC	TIMER expired	Unit switched off after timer expired (→ no error, user information only).
dIFF	Differential Alarm	The differential alarm was triggered The temperature drop in the probe was too quick.
PFAI	Probe failure	The Pt100 probe is broken or damaged The probe was disconnected during operation.
trIA	Triac failure	The heating plate output stage is damaged.
HFAI	Plate Sensor Fail	The hotplate sensor is broken
HOU4	Housing overtemp	The temperature inside the unit has reached its maximum allowed temperature or the internal temperature sensor is damaged
nLI9	Out Of Liquid	The external temperature probe is not inserted into the medium.
COnt	Contact Thermometer Error	The contact thermometer is damaged or was disconnected during operation.
Udt0	Watch Dog Timer	The internal Watch Dog Timer was triggered.
ICEF	Internal Communication Error Front	There was an error in the communication between the printed circuit boards inside the unit
ICE7	Internal Communication Error Motor	There was an error in the communication between the printed circuit boards inside the unit

13.1. Other Messages

LED Display:	Error name:	Description/Reasons :
FAhr	Fahrenheit	°Fahrenheit Mode
OFF	OFF	Unit was switched off by ON/OFF key
4850	RS232 Off	Unit was switched off by RS232/RS485 PC command
HOt	Hotplate is hot	The device is switched off but the hotplate exceeds 60°C.

14. Serial RS485 Interface


The RS485 interface of the Magnetic Hotplate Stirrers M 21, M 22 and M 23 allows remote access of all functions as well as readout and change of all system parameters (e.g. hotplate temperature, motor speed etc.)

Multiple units may be connected parallel on a RS 232 port of a PC with the adapter RS 485 to RS 232 (PN: 61703-0000) and connection cords for connection of different devices via RS 485 (only on demand).

Therefore it is necessary to address the units. Every unit has a so-called slave address between 1 and 255 (default: 1) for parallel operation, each unit needs its own address, to avoid data clashes.

14.1. Pin assignment of the RS485 Interface

The serial port is a RS485 4-wire system with 2 differential inputs (A and B) and 2 differential outputs (Y and Z). The RS485 connector on the back of the unit has the following pin assignment:


Pin 1: RS485 B
 Pin 2: RS485 A
 Pin 3: GND
 Pin 4: reserviert
 Pin 5: reserviert
 Pin 6: +5V (maximal 100mA)
 Pin 7: RS485 Z
 Pin 8: RS485 Y

(RS485 Pin assignment)

	<p>NOTICE</p> <p>Never connect the RJ45 jack with a network jack of a PC, router, switch or something like that. Only connect original accessories to the „RS485 In“ jack! Any other use may result in damage of the device or damage to third party devices!</p>
---	--

14.2. Connection between PC and the device unit


For connecting the unit with a PC a separate RS485 to RS232 adapter (PN: 61703-0000) is necessary. You have to connect the RJ45 plug of this adapter with the „RS485 In“ jack of the unit. The 9-pin Sub-D connector of the adapter is for the connection to the PC and has the following pin assignment:

Pin	Description
2	transmit data TxD
3	receive data RxD
5	ground , GND

For the extension of the cable you can use a standard 9-pin cable (pins not crossed, PN: 30275-0051). When your PC has no RS232 port you can use additionally an USB to RS232 converter (PN: 30244-0001).

Connection between PC and the stirring unit:

NOTICE

Never connect the RJ45 jack with a network jack of a PC, router, switch or something like that. Only connect original accessories to the „RS485 In“ jack! Any other use may result in damage of the device or damage to third party devices!

14.3. Settings of the RS485 Interface

The parameters of the interface (default settings: 9600, 8, N, 1):

Baudrate:	<ul style="list-style-type: none"> • 1200 Baud • 2400 Baud • 4800 Baud • 9600 Baud (default setting)
Databits:	<ul style="list-style-type: none"> • 8 Bit
Parity:	<ul style="list-style-type: none"> • none
Stopbits:	<ul style="list-style-type: none"> • 1
Slave Address:	<ul style="list-style-type: none"> • 1 ... 255 (1 = default setting)

14.4. Addressing of the device unit

To use several units on a bus-system, every unit needs its own address number. This address number is used to send the commands to the unit. No address number must occur twice in the bus system to avoid data collision. The address numbers are possible between 1 and 255. The default setting is 1 for all units. To change the address number you have two possibilities:

Addressing via front panel

You can change the actual RS485 address of the unit easily via the front panel. For changing the address you have to do the following steps:

Press **Power On/Stby** to switch on the device.

Shortly **after** switching the device on **press the Timer Key and hold**.

The display now shows a flashing Γ and the actual RS485 address. Now you can change the address by turning the encoder wheel. Press the encoder wheel to confirm the selected address. The unit will remember the selected address until it is changed again by the user.

Addressing via RS485 command

Alternatively to the addressing via front panel, you have the option to change the slave address with the command WSA:

- Switch the unit on
- Set the new slave address with the RS485 command WSA (example for changing the slave address from 1 to 12: 1,WSA,12<CR>)
- Switch the unit off
- The new address is stored even if the unit is switched on and off

NOTICE

Attention should be paid that the old and the new address is not the same as the slave address of other units in the bus system to avoid data collisions. We recommend that only one unit is connected with the master during the change of an address.

14.5. RS485 Data transfer

Every data transfer is started from the Master (normally the PC). The Master sends a command to the unit (Slave). The slave sends an answer back: a. repeat of the command and the handshake

14.6. Format of a RS485 command

Every command corresponds to the following format:

ADR,CMDCODE,PARAMETERLIST<CR>

Description:

- **ADR:** Slave Address of the desired unit
- **CMDCODE:** Command-code
- **PARAMETERLIST:** 1 to 6 parameters separated by commas
- **CR:** The command string must be terminated by Carriage/Return (ASCII code 13)

14.7. Format of a RS485 Handshake

When a unit gets a command it will answer in this way:

- The complete command is send back
- The unit sends a handshake:

ADR,"HS",RETCODE,PARAMETERLIST<CR>

Description:

- **ADR:** Slave Address of the Controller sending the handshake
- **RETCODE:** Error-code (see table 2 below)
- **PARAMETERLIST:** 0 to 6 parameters (see table 1), each parameter is separated by a comma ","
- **CR:** the handshake as any command, is terminated by ASCII-code 13 (CR)

Return Code	Explanation	Parameterlist
OK	command executed, no error	see table 1
UC	unknown command	None
PA	wrong parameter number (too few or too many parameters specified)	None
NA	command is not allowed in actual operation mode	actual operation mode
PR	at least one parameter is out of range	None
PL	at least one parameter is too long	None
DF	unknown data format	None

NOTICE

The PC must not send the next command until the complete answer is received (repeated command and handshake)!

14.8. RS485 Commands

Overview of all available RS485 Commands:

Command	Function	Parameter Description	Parameter Range	Example	Comment
RTY	Read Type and Version of device	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. name/type of device 2. Version number of software 3. On Off counts 4. Total Minutes Operation time	1 text number number number	1,RTY,1	
PON	Switch on Device	1. Security parameter 1234	1234	1,PON,1234	
OFF	Switch off Device	1. Security parameter 1234	1234	1,OFF,1234	
WON	Set Status (ON/OFF control) of - Motor - Plate	1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	0/1 0/1	1,WON,1,0	
RON	Read status of: - Motor - Plate	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	1 0/1 0/1	1,RON,1	
RAC	Read actual	1. Dummy parameter to initiate transfer -> Controller sends in handshake 1. actual motor speed 2. actual plate temperature in °C 3. actual probe temperature in °C (if connected, if not 'x' as return value) 4. Parameter reserved 5. last off-Condition	1 0..1600 0..MAX 0..250, x x 100...142* see: table RS485 Off condition numbers	1,RAC,1	Motor disabled in H30/30D. Device will send 'x' as parameter.
WSE	Write actual set values of: - Motor - Hotplate (plate) - external PT100 sensor (probe)	1. Setpoint of motor speed in Rpm 2. Setpoint of hotplate temperature in °C 3. Setpoint of probe temperature in °C	0..1600 0..MAX platetemp 0..250	1,WSE,800,340,60	MAX depends on plate-type (Ceramic-500°C, Aluminium 380°C, Stainless steel 350°C, Heatable Flask Stirrer 450°C) Motor disabled in H30/30D. Device will ignore the value.
RSE	Read actual set values of: - Motor - Hotplate (plate) - external PT100 sensor (probe)	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Setpoint of motor speed in Rpm 2. Setpoint of hotplate temperature in °C 3. Setpoint of probe temperature in °C	1 0..1600 0.. MAX platetemp 0..250	1,RSE, 1	Motor disabled in H30/30D. Device will send 'x' as parameter.
WTR	Set / Write: -Timer	1. Timer value in seconds. Set this value to 0 to disable the timer.	0.. 59940	1,WTR,600,450,160	Ramp disabled in M 21, M 22, M 23 KM16.4, KM16.7 and

	- Ramp - Safety temperature	2. setting of the ramp in °C/h (a value of 450 disables the ramp) 3. Safety temperature in °C	1.. 450 20.. MAX platetemp+25°C		H30/30D. Device will ignore the value.
RTR	Read settings of: - Timer - Ramp - Safety temperature	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. actual setting of the timer in seconds. If the timervalue is 0, the timer is disabled. 2. setting of the ramp in °C/h (a value of 450 signals that the ramp is disabled) 3. Safety temperature in °C	1 0.. 59940 1..450, x 20.. MAX platetemp+25°C	1,RTR,1	Ramp disabled in M 21, M 22, M 23 KM16.4, KM16.7 and H30/30D. Device will send 'x' as parameter.
WVO	Write/Set volume	1. Volume in ml	100..10000	1,WVO,1000	
RVO	Read volume	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. programmed volume in ml	1 100..10000	1,RVO,1	
WSM	Activate serial mode (→ enable/disable setpoint control via frontpanel)	1. Set this value to 0 to disable the motor setpoint control via RS232 commands, set this value to 1 to enable motor setpoint control via RS232 commands and disable the frontpanel control. (only used for units with motor)	0 / 1		
RTU	Read Temperature Units	1. type of units	0=Celsius / 1=Fahrenheit	1,RTU,1	
RSS	Read system status	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Device On/Standby 0 → Standby 1 → On 2. Reserved, device will send 'x' as parameter	0/1 x	1,RSS,1	
WSA	Set RS485 slave address	1. New slave address of controller	0...255	1,WSA,3	
WBD	Set RS485 Baudrate	1. Baudrate	0=1200baud / 1=2400baud / 2=4800baud / 3=9600baud	1,WBD,2	After setting the new baudrate you have to switch off and on the device to store the new value in EEPROM

*Table RS485 Off condition numbers

Off condition Number	Off Condition	Description
100	NOOFFCONDITION	Undefined Off Condition
101	OFF_USER	Unit was switched off by ON/OFF key
102	OFF_REMOTE	Unit was switched off by RS232/RS485 PC command
103	TIMER_EXPIRED	Unit switched off after timer expired
106	PROBE_DIFFERENTIALALARM	The differential alarm was triggered The temperature drop in the probe was too quick
107	PROBE_OUTOFLIQUID	The external temperature probe is not inserted into the medium
108	PROBE_SAFETY	The measured probe plate temperature was higher than the set safety temperature for the Pt100 probe (Safety Temp.)
113	PROBE_BROKEA	The Pt100 probe is broken or damaged The probe was disconnected during operation.
117	CONTACTTHERMOMETER_BROKE	The contact thermometer is damaged or was disconnected during operation.
118	PLATE_OVERTEMP	The hotplate temperature exceeded its maximum value The hotplate temperature sensor may be broken
120	PLATE_SAFETY	The measured hotplate temperature was higher than the set safety temperature for the hotplate (Safety Temp.)
125	PLATE_BROKEA	The hotplate sensor is broken
130	PLATE_SHORTEDHOTPLATE	The temperature rise in the hotplate was too fast. The heating plate output stage is damaged.
134	DEVICE_INTERNALCOMMUNICATIONERROR_FRONT	There was an error in the communication between the printed circuit boards inside the unit
135	DEVICE_INTERNALCOMMUNICATIONERROR_MOTOR	There was an error in the communication between the printed circuit boards inside the unit
139	DEVICE_INTERNALTEMPERROR	The temperature inside the unit has reached its maximum allowed temperature or the internal temperature sensor is damaged
140	DEVICE_POWERERROR	The device was not switched off with the power button. Disconnecting the power cable or a mains power failure may be the source for error.
141	DEVICE_WDTERORR	The internal Watch Dog Timer was triggered

14.9. RS485 Examples of use

Example 1: Setting of values

Programming of a M 21 unit with the following parameters:

- Hotplate temperature: 300°C (corresponds to the plate limit)
- Probe temperature: 50°C
- Motor speed: 500rpm

The unit has the standard Slave Address 1, a PT100 sensor is connected:

Following commands must be sent to the unit, the answers of the M 21 are in *italic types*, they must not be sent to the unit:

- **1,PON,1234<CR>** ;Switch the unit on
 - *1,PON,1234<CR>* ;Command is sent back repeatedly from the device
 - *1,HS,OK<CR>* ;Handshake sent from the unit, command accepted
- **1,WSM,1<CR>** ;deactivating of the front panel only controlling via RS485 (optional)
 - *1,WSM,1<CR>* ;Command is sent back repeatedly from the device
 - *1,HS,OK<CR>* ;Handshake sent from the unit, command accepted
- **1,WSE,500,300,50<CR>** ;set values (motor speed, hotplate and Probe temperature)
 - *1,WSE,500,300,50<CR>* ;Command is sent back repeatedly from the device

- **1,HS,OK<CR>** ;Handshake sent from the unit, command accepted
- **1,WON,1,1<CR>** ;Switch on motor and hotplate
 - **1,WON,1,1<CR>** ;Command is sent back repeatedly from the device
 - **1,HS,OK<CR>** ;Handshake sent from the unit, command accepted

Example 2: Reading of actual values

Reading of the actual values of a M 21 unit:

- Temperature units
- Hotplate temperature
- Probe temperature
- Motor speed

The unit has the standard Slave Address 1, a PT100 sensor is connected:

Following commands must be send to the unit, the answers of the M 21 are in *italic types*, they must not be send to the unit:

- **1,PON,1234<CR>** ;Switch the unit on
 - **1,PON,1234<CR>** ;Command is sent back repeatedly from the device
 - **1,HS,OK<CR>** ;Handshake sent from the unit, command accepted
- **1,RTU,1<CR>** ;Reading temperature units
 - **1,RTU,1<CR>** ;Command is sent back repeatedly from the device
 - **1,HS,OK,0<CR>** ;Handshake sent from the unit, command accepted
(units = Celsius)
- **1,RAC,1<CR>** ;Reading actual Values
 - **1,RAC,1<CR>** ;Command is sent back repeatedly from the device
 - **1,HS,OK,480,180,50,x,4<CR>** ;Handshake sent from the unit, command accepted
(motor speed = 480rpm, hotplate temp = 180°C, probe temp = 50°C)

The actual settings are shown in the parameter lists:

- The temperatures are shown in Celsius (command RTU, parameter 1 = 0)
- The actual motor speed is 480rpm (command RAC, parameter 1 = 480)
- The actual hotplate temperature is 180°C (command RAC, parameter 2 = 180)
- The actual probe temperature is 50°C (command RAC, Parameter 3 = 50)

15. Maintenance and Cleaning

▲WARNING

In case of malfunction do not attempt to repair the unit. There are no user-serviceable parts in this instrument.

The outer casting is mostly 304 grade stainless steel, the hotplate consist of stainless steel (M 21), anodized aluminium (M 22) or Ceran (M 23). Together with a chemically resistant splash-proof membrane key-pad and therefore easily cleaned with warm water and any proprietary liquid laboratory detergent. Do not use steel wool or any similar plastic wool sponge or any aggressive cleaning agents to clean the unit.

Crust, limestone and water stains may be cleaned with the usual household detergents. Please remove remaining detergent thoroughly with a damp cloth or sponge, since some of the detergents react caustic with higher temperatures. If soils are already encrusted use a razorblade scratcher. Plastics, kitchen foil as well as stuff containing sugar must be removed immediately.

Mechanical supports are: Razorblade scratchers, non-scratching sponges and non-scratching steel wool.

The M 21, M 22 and M 23 should only be opened and repaired by authorized service personnel. Any work on the electronics in the unit should only be carried out by knowledgeable, trained personnel. Any attempt by the user to repair the unit will immediately render the guarantee null and void. Please contact your local distributor in the event of a problem.

16. Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's election, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement unit will be warranted for the remainder of the original twenty-four (24) month period of ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein. The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damage of any kind resulting from use or possession of the product.

▲WARNING

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

17. Dismantling, Transport and Storage

17.1. Dismantling

1. Switch off the unit by switching the rocker switch (1) to the „OFF“-Position.
2. Set speed to “0” position by turning the speed control knob to “0”.
3. Disconnect the unit from the mains.
4. Now you may remove the instrument from the working area.

17.2. Transport and Storage

Prior to transport:

Switch the instrument off and unplug the power supply.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, please sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

WEEE No. DE 54398924

18. Technical Data

Type M 21, M 22, M 23	Specifications
Mains	115/230 Volt; 50/60 Hz (See rear panel)
Display	4 digit 7Segment-LED display
Motor	– split-pole, 25 Watts – 60 - 1600 min ⁻¹ – Softstart, microprocessor controlled

Temperature control	<ul style="list-style-type: none"> – microprocessor controlled (PID Control) – Pt100 probe adjustable in 1°C Steps – Hotplate adjustable in 1°C steps up to 350°C (M 21), 380°C (M 22), 500°C (M 23) – programmable probe temperature (Pt100) – programmable hotplate temperature – programmable safety temperature
Safety system	<p>detects and protects against the following hazardous situations:</p> <ul style="list-style-type: none"> -hotplate failure (hotplate exceeds max. allowed temperature) -probe failure (disconnection or break of Pt100 probe) -detection of out-of-liquid condition, e.g. stirring vessel breaks -differential safety cut-out (external probe senses an extreme drop in temperature of liquid being stirred/heated) -triac failure
Switch-off time delay (Timer)	<p>programmable switch-off time delay</p> <p>min. delay: 1 minute</p> <p>max. delay: 999 minutes</p>
Temperature sensors	<p>Hotplate: Thermocouple</p> <p>External: Platinum Resistance Sensor (Pt100)</p>
Setpoint resolution hot plate	1°C / 1K (1°F in Fahrenheit Mode)
Setpoint resolution Pt100 Probe	1°C / 1K (1°F in Fahrenheit Mode)
Motor speed resolution	10 min ⁻¹
Typical temperature stability with Pt100 use*	± 0.4 °C*
Temperature range heating medium	Up to 250°C
Recommended stirring bar	<p>Length: 55 mm (30 – 70 mm)</p> <p>Diameter: 10 mm (5 – 12 mm)</p>
Permissible Ambient Temperature	5 - 40 °C
Permissible Humidity	80 % RH
Safety class acc to DIN 40050	IP32
Case dimensions W x D x H	145 x 210 x 110 mm
Weight	2,6kg

Type M 21	
Hotplate material	Stainless steel
Hotplate dimensions	∅ 135mm round
Heating power	500 Watts

Type M 22	
Hotplate material	Anodised aluminium
Hotplate dimensions	∅ 140mm round
Heating power	500 Watts

Type M 23	
Hotplate material	Ceran
Hotplate dimensions	135mm square
Heating power	600 Watts

* 2L H₂O at 60°C Setpoint, Ambient: 23°C. Maximum variance during 60 min test.

19. Repair Return Form

<p>Important:</p> <ul style="list-style-type: none"> - Calibration Service is automatically performed for instruments that require adjustments - We contact you by sending a repair estimate. - Repairs will only be carried out <u>after</u> your release.
--

CONTACT/USER INFORMATION	
Contact:	Phone No.:
Fax No.:	Email:
Billing:	Shipping:
Company	Company
Address	Address

INSTRUMENT INFORMATION	
Model	Serial No.
Please describe all problems/malfunctions	

Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to Ingenieurbüro CAT, M. Zipperer GmbH and an MSDS for hazardous material must be included with the instrument.

RETURN SHIPPING			
UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>

**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

PACKAGING INSTRUCTIONS TO RETURN AN INSTRUMENT FOR REPAIR

- Remove all accessories (e.g. homogeniser tools, stirring paddles, stands, clamps, cables) from the instrument
- Clean excess testing material off the instrument/accessory
- Include MSDS sheets for all hazardous materials used with this instrument
- Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- DO NOT send accessories unless there is a problem with them. When sending back, wrap each item to avoid contact with the instrument.
- Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.
Example: "Handle with care" or "Fragile- Delicate Instrument" and send to:

Ingenieurbüro CAT
 M. Zipperer GmbH
 Wettelbrunner Str. 6
 D-79282 Ballrechten-Dottingen

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Command	Function	Parameter Description	Parameter Range	Example	Comment
RTY	Read Type and Version of device	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. name/type of device 2. Version number of software 3. On Off counts 4. Total Minutes Operation time	1 text number number number	1,RTY,1	
PON	Switch on Device	1. Security parameter 1234	1234	1,PON,1234	
OFF	Switch off Device	1. Security parameter 1234	1234	1,OFF,1234	
WON	Set Status (ON/OFF control) of - Motor - Plate	1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	0/1 0/1	1,WON,1,0	
RON	Read status of: - Motor - Plate	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	1 0/1 0/1	1,RON,1	
RAC	Read actual	1. Dummy parameter to initiate transfer -> Controller sends in handshake 1. actual motor speed 2. actual plate temperature in °C 3. actual probe temperature in °C (if connected, if not 'x' as return value) 4. Parameter reserved 5. last off-Condition	1 0..1600 0..MAX 0..250, x x 100...142* see: table RS485 Offcondition numbers	1,RAC,1	Motor disabled in H30/30D. Device will send 'x' as parameter.
WSE	Write actual set values of: - Motor - Hotplate (plate) - external PT100 sensor (probe)	1. Setpoint of motorspeed in Rpm 2. Setpoint of hotplate temperature in °C 3. Setpoint of probe temperature in °C	0..1600 0..MAX platetemp 0..250	1,WSE,800,340,60	MAX depends on plate-type (Ceramic- 500°C, Aluminium 380°C, Stainless steel 350°C, Heatable Flask Stirrer 450°C) Motor disabled in H30/30D. Device will ignore the value.
RSE	Read actual set values of: - Motor - Hotplate (plate) - external PT100 sensor (probe)	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Setpoint of motorspeed in Rpm 2. Setpoint of hotplate temperature in °C 3. Setpoint of probe temperature in °C	1 0..1600 0.. MAX platetemp 0..250	1,RSE, 1	Motor disabled in H30/30D. Device will send 'x' as parameter.

WTR	Set / Write: -Timer - Ramp - Safety temperature	1. Timer value in seconds. Set this value to 0 to disable the timer. 2. setting of the ramp in °C/h (a value of 450 disables the ramp) 3. Safety temperature in °C	0.. 59940 1.. 450 20.. MAX platetemp+25°C	1,WTR,600,450,160	Ramp disabled in M 21, M 22, M 23 KM16.4, KM16.7 and H30/30D. Device will ignore the value.
RTR	Read settings of: - Timer - Ramp - Safety temperature	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. actual setting of the timer in seconds. If the timervalue is 0, the timer is disabled. 2. setting of the ramp in °C/h (a value of 450 signals that the ramp is disabled) 3. Safety temperature in °C	1 0.. 59940 1..450, x 20.. MAX platetemp+25°C	1,RTR,1	Ramp disabled in M 21, M 22, M 23 KM16.4, KM16.7 and H30/30D. Device will send 'x' as parameter.
WVO	Write/Set volume	1. Volume in ml	100..10000	1,WVO,1000	
RVO	Read volume	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. programmed volume in ml	1 100..10000	1,RVO,1	
WSM	Activate serial mode (-> enable/disable setpoint control via frontpanel)	1. Set this value to 0 to disable the motor setpoint control via RS232 comands, set this value to 1 to enable motor setpoint control via RS232 commands and disable the frontpanel control. (only used for units with motor)	0 / 1		
RTU	Read Temperature Units	1. type of units	0=Celsius / 1=Fahrenheit	1,RTU,1	
RSS	Read system status	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Device On/Standby 0 → Standby 1 → On 2. Reserved, device will send 'x' as parameter	0/1 x	1,RSS,1	
WSA	Set RS485 slave address	1. New slave address of controller	0...255	1,WSA,3	
WBD	Set RS485 Baudrate	1. Baudrate	0=1200baud / 1=2400baud / 2=4800baud / 3=9600baud	1,WBD,2	After setting the new baudrate you have to switch off and on the device to store the new value in EEPROM

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