



Ingenieurbüro CAT
M. Zipperer GmbH

Instruction Manual

Gebrauchsanleitung

Homogenising System
Dispergiersystem

UNIDRIVE X 1000 X 1000D



Content

| | | |
|-------|---|----|
| 1 | User Instructions | 5 |
| 1.1 | Important Instructions for your safety..... | 5 |
| 1.2 | Danger symbols in this operating manual..... | 5 |
| 1.2.1 | Danger symbols:..... | 5 |
| 1.2.2 | Danger levels..... | 6 |
| 1.2.3 | Warning signs on the device | 6 |
| 2 | General safety warnings and instructions | 7 |
| 3 | General Information | 9 |
| 4 | Intended Use | 9 |
| 5 | Scope of delivery..... | 10 |
| 5.1 | Homogenising Tools and Accessories (to be ordered separately):..... | 10 |
| 6 | Description Drive Unit UNIDRIVE X 1000/X 1000D..... | 10 |
| 7 | Setting up and Starting up the Instrument..... | 12 |
| 7.1 | Unpacking the Instrument | 12 |
| 7.2 | Setting up the Instrument | 12 |
| 7.3 | Electrical Connection..... | 12 |
| 7.4 | Assembly of Support Rod | 13 |
| 7.5 | Attaching the Drive Unit to a Stand | 14 |
| 7.6 | Description of a Homogenising Tool | 14 |
| 7.7 | Inserting the Homogenizing Tool | 15 |
| 8 | Operation UNIDRIVE X 1000/X 1000D | 16 |
| 8.1 | Connect the Unit to the Mains..... | 17 |
| 8.2 | Description of the Speed Scale X 1000 and X 1000D: | 17 |
| 8.3 | Exchanging the Homogenizing Shafts | 17 |
| 8.4 | Exchanging the Generator:..... | 17 |
| 8.5 | Limits on continuous operation of the X 1000 drive unit | 18 |
| 8.6 | Overheating protection for the X 1000 drive unit | 18 |
| 9 | Maintenance and Cleaning | 19 |
| 9.1 | Cleaning the Drive Unit UNIDRIVE X 1000/X 1000D | 19 |
| 9.2 | Cleaning the Homogenising Tool..... | 19 |
| 9.3 | Rough Cleaning of Flow-Through-Chamber and Tubing | 19 |

| | | |
|------|--|----|
| 9.4 | Cleaning the Inside of the Flow-Through-Chamber | 19 |
| 9.5 | Maintenance of drive units UNIDRIVE X 1000/X 1000D | 20 |
| 9.6 | Maintenance homogenizing tools..... | 20 |
| 10 | Dismantling, Transport and Storage | 20 |
| 10.1 | Dismantling | 20 |
| 10.2 | Transport and Storage | 20 |
| 11 | Disposal | 21 |
| 12 | Warranty and Liability..... | 21 |
| 13 | Technical Data..... | 22 |
| 13.1 | Technical Data UNIDRIVE X 1000/X 1000D | 22 |
| 14 | Repairs..... | 22 |
| 14.1 | Repair Return Form | 23 |

Inhaltsverzeichnis

| | | |
|-------|---|----|
| 1 | Allgemeine Information | 24 |
| 2 | Auspicken des Gerätes..... | 24 |
| 2.1 | Lieferumfang und Zubehör..... | 24 |
| 2.2 | Dispergierwerkzeuge und Zubehör (separat zu bestellen):..... | 25 |
| 3 | Sicherheitshinweise | 27 |
| 3.1 | Erläuterung der Sicherheitshinweise in der Gebrauchsanleitung | 27 |
| 3.2 | Erläuterung der Sicherheitshinweise auf dem Gerät | 27 |
| 3.3 | Allgemeine Sicherheitshinweise..... | 27 |
| 4 | Aufbau und Betrieb | 28 |
| 4.1 | Zweckbestimmung | 28 |
| 4.2 | Beschreibung des Antriebs UNIDRIVE X 1000/X 1000D..... | 29 |
| 4.3 | Aufbau | 30 |
| 4.3.1 | Montage der Haltestange | 30 |
| 4.3.2 | Befestigen des Antriebs am Stativ | 31 |
| 4.4 | Beschreibung eines Dispergierwerkzeugs..... | 31 |
| 4.4.1 | Anbringen des Dispergierwerkzeugs | 32 |
| 4.5 | Betrieb UNIDRIVE X 1000/X 1000D | 33 |

| | | |
|-------|---|----|
| 4.5.1 | Anschließen des Gerätes ans Stromnetz | 34 |
| 4.5.2 | Beschreibung der Stellradzahlen des X 1000 und X 1000D: | 34 |
| 4.6 | Hinweise zum Dauerbetrieb der X 1000-Dispergierantriebe..... | 34 |
| 4.7 | Überhitzungsschutz der X 1000-Dispergierantriebe | 34 |
| 4.8 | Austausch der Dispergierwerkzeuge | 35 |
| 4.8.1 | Austausch des Generators | 35 |
| 5 | Reinigung und Wartung | 36 |
| 5.1 | Reinigung des Antriebs UNIDRIVE X 1000/X 1000D..... | 36 |
| 5.2 | Reinigen des Dispergierwerkzeugs..... | 36 |
| 5.3 | Vorreinigung der Durchflusskammer und den Schläuchen | 36 |
| 5.4 | Reinigung des Kammerinnenraums | 36 |
| 5.5 | Wartung der Dispergierwerkzeuge | 37 |
| 6 | Abbau, Transport und Lagerung | 37 |
| 6.1 | Abbau | 37 |
| 6.2 | Transport und Lagerung | 37 |
| 7 | Entsorgung | 38 |
| 8 | Garantie und Haftungsausschluss..... | 38 |
| 9 | Reparaturen | 38 |
| 10 | Technische Daten | 39 |
| 10.1 | Rücksendeformular..... | 40 |
| 11 | Declaration of conformity / Konformitätserklärung..... | 41 |

1 User Instructions

1.1 Important Instructions for your safety



- Every user must read and understand this manual completely before use. Failure to do so can result in serious injury or death.
- Comply with all safety and accident-prevention regulations applicable to laboratory work.
- Follow general instructions for hazard prevention and general safety instructions, e.g. wear protection clothing, eye protection and gloves.
- This operating manual is part of the product. Thus, it must always be easily accessible.
- This instruction sheet does not purport to address all of the safety problems which might result from the use of this device, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult their authorized safety advisors and establish appropriate health and safety practices and then determine the application of regulatory limitations prior to use.
- Enclose this operating manual when transferring the device to another place.
- If this manual is lost, please request another one. Please contact your dealer or

Ingenieurbüro CAT
M. Zipperer GmbH
Wettelbrunner Str. 6
D-79282 Ballrechten-Dottingen
Tel.: ++49-(0)7634-5056-800
Fax: ++49-(0)7634-5056-801
www.cat-ing.de
info@cat-ing.de

1.2 Danger symbols in this operating manual

The safety instructions in this manual appear with the following danger symbols and danger levels:

1.2.1 Danger symbols:

| | | | |
|--|--------------|---|------------------|
| A warning symbol consisting of an exclamation mark inside an equilateral triangle with a black border. | Hazard point | A warning symbol consisting of a lightning bolt inside an equilateral triangle with a black border. | Electrical shock |
|--|--------------|---|------------------|

| | | | |
|--|-----------------|--|-----------------|
| | Risk of fire | | Explosion |
| | Crushing | | Hot surface |
| | Bio hazard | | Magnetic field |
| | Chemical hazard | | Material damage |

1.2.2 Danger levels

| | |
|----------------|--|
| DANGER | Will lead to severe injuries or death |
| WARNING | May lead to severe injuries or death |
| CAUTION | May lead to light to moderate injuries |
| NOTICE | May lead to material damage |

1.2.3 Warning signs on the device

| | |
|--|---|
| | WARNING These symbols indicate that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user. |
| | WARNING These symbols indicate that it is imperative to wear eye and ear protection to avoid damage to health |

2 General safety warnings and instructions

| | |
|---|---|
|  | !DANGER Risk of explosion. <ul style="list-style-type: none">• Do not operate the device in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof.• Do not use this device for processing any substances which could generate an explosive atmosphere.• Do not use this device to process any explosive or highly reactive substances.• Do not use flammable or explosive substances near the instrument. |
|  | !DANGER Electric shock as a result of penetration of liquid. <ul style="list-style-type: none">• Do not allow any liquids to penetrate the inside of the housing• Switch off the device and disconnect the power plug before starting cleaning or disinfection work. The On/Off Switch on the device does not disconnect the device from the power source.• Use only closed tubes.• Only plug the device back in if it is completely dry, both inside and outside. |
|  | Risk from incorrect supply voltage <ul style="list-style-type: none">• Only connect the device to voltage sources which correspond to the electrical requirements on the type label. |
|  | !WARNING Electric shock due to damage to device or mains cable <ul style="list-style-type: none">• Only connect the device to the mains supply if the device and the mains cable are undamaged• Only use devices that have been properly installed or repaired.• In case of danger, disconnect the device from the mains supply by pulling the power plug from the mains socket or by using the isolating device intended for this purpose (e.g. emergency stop switch) |
|  | !WARNING Lethal voltage inside the device <ul style="list-style-type: none">• Do not open the device.• Ensure that the housing is always closed and undamaged so that no parts inside the housing can be contacted by accident.• The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.• Do not allow any liquids to penetrate the inside of the housing.• Repairs are only to be carried out by trained service technicians. |



⚠WARNING **Damage to health due to corrosive or aggressive chemicals**

- Observe all markings on the reagent bottles.
- Always check the instrument that all push-ons and threaded connections are firmly in place before beginning operation.
- Use proper connecting vessels, protective clothing and gloves.
- Avoid splashes
- When dispensing, maintain a physical distance between the instrument and the body.
- Dangerous and fuming chemicals must be dispensed in a fume hood.
- Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction.
- Always use the instrument in such a manner that neither the operator, nor any other person is endangered.



⚠WARNING **Damages to health due to infectious liquids and pathogenic germs.**

- When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets and the manufacturer's application notes.
- Wear personal protective equipment
- For comprehensive regulations about handling germs or biological material of the risk group II or higher, please refer to the "Laboratory Biosafety Manual" in its respectively current valid version from the World Health Organisation



⚠WARNING **Damage to health due to contaminated device and accessories**

In the following cases, sample material can be released:

- improperly sealed tubes
 - unstable tubes
 - high vapour pressure of the content so that the seal of the tubes can spring open
 - damaged sealings
 - smashed glass tubes
- Only mix in closed tubes
 - Observe the nationally prescribed safety environment when working with hazardous, toxic and pathogenic samples. Pay particular attention to personal protective equipment (gloves, clothing, goggles, etc.), extraction, and the safety class of the lab.
 - Decontaminate the device and the accessories before storage and shipping.



⚠WARNING **Risk of fire**

- Do not use this device to process any highly flammable liquids

| | |
|---|---|
|  | <p>⚠ CAUTION Poor safety due to inadequate fixing of the unit</p> <ul style="list-style-type: none"> • Ensure that the unit is firmly attached to a solid stand. |
|  | <p>⚠ CAUTION Poor safety due to incorrect accessories and spare parts.</p> <p>The use of accessories and spare parts other than recommended by Ingenieurbüro CAT, M. Zipperer GmbH may impair the safety, function and precision of the device. Ingenieurbüro CAT, M. Zipperer GmbH cannot be held liable or accept any liability for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.</p> <ul style="list-style-type: none"> • Only use accessories and spare parts recommended by Ingenieurbüro CAT, M. Zipperer GmbH |
|  | <p>⚠ CAUTION Crush hazard due to moving parts</p> <ul style="list-style-type: none"> • Do not replace any consumables as long as the device is running. • Do not open the coverage as long as the device is running |

3 General Information

The homogeniser Unidrive X 1000/Unidrive X 1000D is designed in accordance with Safety Class 1 and built and tested in accordance with DIN EN 61010.

According to these regulations the unit is designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument the user should follow the instructions and safety guidelines in this manual.

4 Intended Use

The UNIDRIVE X 1000/X 1000D is a homogenizing system which in connection with a homogenizing tool produces emulsions and dispersions. For safe operation the unit has to be fixed to a stand.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.


| | |
|---|--|
|  | <p>⚠ WARNING It is the user's task to find out whether the device is suitable for the application. If in doubt clarify this with your dealer or directly with the manufacturer.</p> |
|  | <p>⚠ WARNING Please comply with all safety and accident-prevention regulations applicable to laboratory work.</p> |

5 Scope of delivery

Please check that the package contains the following:

| Type | Description | Part No. |
|---------|-------------------------------|---------------|
| X 1000 | Homogenizing Drive Unit, 230V | 60402-0050 or |
| X 1000 | Homogenizing Drive Unit, 115V | 60402-0051 or |
| X 1000D | Homogenizing Drive Unit, 230V | 60406-0050 or |
| X 1000D | Homogenizing Drive Unit, 115V | 60406-0051 or |
| | 1 Support Rod | |
| | 1 Instruction Manual | |

5.1 Homogenising Tools and Accessories (to be ordered separately):

| T 6 | T 6 long | T 10 | T 17 | T 20 | G 20 |
|---|---|--|---|---|---|
|  |  |  |  |  | |
| 2,0 - 50 ml | 10,0 - 250 ml | 250 - 1500 ml | 250 - 1500 ml | 250 - 5000 ml | 250 - 5000 ml |
| 55 mm | 120 mm | 150 mm | 150 mm | 200 mm | 215 mm |
|  |  |  |  |  | |
| 60420-000F | 60410-000F | 60421-000F | 60426-000F | 60422-000F | 60423-000F |
|  |  |  |  |  | |
| 60420-000N | 60410-000N | 60421-000N | 60426-000N | 60422-000N | 60423-000N |
|  |  |  |  |  | |
| 60420-000V | 60410-000V | 60421-000V | 60426-000V | 60422-000V | 60423-000V |
|  |  |  |  |  | |
| 60420-000M | 60410-000M | 60421-000M | 60426-000M | | 60423-000M |

T = Shaft with PTFE bearing. This shaft type is used for standard application. (For liquid media like water, oil, etc.)

The bearing must be cooled by the media, *which penetrates into the bearing*.


A dry run is not permitted.

G = Shaft with slide seal ring. This shaft type is excellent for applications with aggressive respectively abrasive media and when pressure is applied. *Also suitable for use with solids.*

The bearing is sealed and suitable for dry running

| Tool | Description | Part No. |
|--------------------------|--|------------|
| Pedestal stand incl. rod | platform overall dimensions: 330x200 mm, rod diameter: 16mm, length:690mm | 60491-0000 |
| Cross over clamp | Casted alloy clamp to attach the unit to a stand | 60492-0000 |
| Cross over clamp | Steel clamp to attach the unit to a stand | 60492-0050 |
| Universal key | Universal rotor wrench for Shafts 6 - 30 mm | 60470-0000 |
| Socket wrench | for all shafts | 60471-0000 |
| Rotor wrench | for shaft T 10 | 60472-0000 |

6 Description Drive Unit UNIDRIVE X 1000/X 1000D



7 Setting up and Starting up the Instrument

7.1 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

If the instrument is not damaged and all parts are complete you may start to operate the device after reading the instruction manual.

7.2 Setting up the Instrument


| | | |
|---|----------------|--|
|  | DANGER | The device may not be operated in explosion-prone areas. |
|  | WARNING | The device is not to be used without supervision. |
|  | WARNING | When connecting the instrument to an AC power outlet ensure that your local supply voltage matches the indication on the instrument. |

7.3 Electrical Connection

WARNING

- On delivery, the units are supplied with a Euro plug EN 50075:1990, 2.5 A, 250 V. For North America, with a U.S. standard plug (NEMA 1, Type A) 10 A, 125 V.
- When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.
- When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

7.4 Assembly of Support Rod




- Screw the support rod (8) into the threaded bore (7) at the back of the unit.
- Insert a screw driver into the bore at the end of the support rod (8) and firmly tighten the support rod by turning the screw driver to the right.

7.5 Attaching the Drive Unit to a Stand


**WARNING**

To ensure safe operation the drive units are mounted to a pedestal stand using a cross over clamp. Vibrations may loosen the holding screws. Please check periodically whether all holding screws are securely tight



- Attach cross over clamp (9) to a stand rod.
- Now attach the drive unit to the stand by inserting the support rod (8) into the free opening of the cross over clamp (9).
- Tighten screws (10) of the cross over clamp (9).

7.6 Description of a Homogenising Tool




7.7 Inserting the Homogenizing Tool



⚠️WARNING

The generators have sharp edges. Handle with care. Risk of injuries!

The homogenizing tool comes always assembled. The connecting adapter (drive pin (11)) is located on the top of the Shaft and is used to attach the Shaft to the drive motor. When inserting the homogenising tool, make sure that the coupling fits well. This can be easily detected by the "click" that is heard when the coupling Shaft has reached the end of the coupling slot of the motor. Securely fasten the holding screw (6).



NOTICE

Never run homogenising tools dry as the bearings and gaskets will be damaged if the generator is not cooled by the medium. The distance between the homogenising tool and the vessel bottom should not be less than 10mm. The filling of the medium may not be less than about 55 mm. To avoid a vortex it is necessary to insert the dispersing shaft out of the centre of the vessel. If these conditions are met the unit is ready for operation.

8 Operation UNIDRIVE X 1000/X 1000D

**⚠WARNING**

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

Connect the unit to the mains.

The drive unit is turned on and off by means of the integral switch (1) at the motor housing. With the integral switch (1) the speed can be steplessly adjusted.

**⚠WARNING**

The integral switch (1) is only single-pole and is thus not suitable for safely isolating the unit from the mains supply!

**⚠WARNING**

Before starting the unit make sure that

- the speed is set to the lowest speed
- the generator of the homogenising tool is inserted in the fluid. To avoid a vortex, it is necessary to insert the dispersing shaft out of the centre of the vessel.
- the vessel is protected against rotating, shaking or moving
- make sure that the air ventilation slots are always open, as this is necessary for cooling purposes

**⚠WARNING**

Start the unit always at the lowest speed and then gradually increase the speed to the desired level.

NOTICE

Never exceed the maximum operation time of 3 minutes in one go. Allow motor and shaft to cool down. Disregard cause serious damage to drive unit and shaft.

NOTICE

The UNIDRIVE X 1000/X 1000D has an internal thermal overload protection which switches off the instrument to protect the motor.


Before restart the unit wait for minimum 10 min. to let the motor cool down.

8.1 Connect the Unit to the Mains.

The drive unit is turned on and off by means of the on/off switch (3) at the front of the motor housing.

- position I: motor runs
- position O: motor off

Set the speed with the speed adjustment wheel (2).



Before starting the unit make sure that

- the speed is set to the lowest speed
- the Generator (12) of the homogenising tool is inserted in the fluid. To avoid a vortex it is necessary to insert the dispersing Shaft out of the centre of the vessel.
- the vessel is protected against rotating, shaking or moving
- make sure that the air ventilation slots (4) are always open, as this is necessary for cooling purposes

Start the unit always at the lowest speed and then gradually increase the speed to the desired level.

8.2 Description of the Speed Scale X 1000 and X 1000D:


| | | |
|---|---|-----------|
| 1 | = | 4000 rpm |
| 2 | = | 8500 rpm |
| 3 | = | 15000 rpm |
| 4 | = | 20000 rpm |
| 5 | = | 25000 rpm |
| 6 | = | 33000 rpm |

The set speed of UNIDRIVE X 1000D is shown the speed display (5).

8.3 Exchanging the Homogenizing Shafts

- Prior to exchange the homogenizing shaft disconnect the unit from the mains supply by unplugging the unit.
- Loosen the holding screw (6) and firmly pull out the Shaft by hand. It will disconnect easily.
- When inserting another Shaft, make sure that the coupling fits well. This can be easily detected by the "click" that is heard when the coupling Shaft has reached the end of the coupling slot of the motor.
- Securely fasten the holding screw (6).

8.4 Exchanging the Generator:



To exchange the Generator (rotor and stator), always take the Shaft out of the motor drive unit as described under 4.6.

For undoing the rotor, use the two special keys (to be ordered separately). The socket wrench is used to counter-hold the Shaft, and the rotor wrench to undo the rotor.

To exchange the stator, remember that the fastening thread is **counter-clockwise**.

You will find further detailed descriptions in the instruction manual which comes with each homogenising tool.

NOTICE

Always unscrew the rotor first before unscrewing the stator.

8.5 Limits on continuous operation of the X 1000 drive unit

In combination with the shafts T6, T10, T17 and T/G20 the X 1000 drive unit can be used in continuous operation.

8.6 Overheating protection for the X 1000 drive unit

The X 1000 drive unit features a safety electronic (temperature-dependent overload protection), which prevent the unit from overheating by enabling a cooling mode as soon as the drive reaches the critical temperature.

In that cooling mode the constant-electronic is deactivated and the drive unit is running at a very low speed. After a cooling time of 60 seconds, the machine is ready for operation. Therefore, switch off and on the device to activate the constant-electronic. Before switching on the device, ensure that the lowest rpm-level is set to prevent the media from splashing out of the container. The speed can then be raised step-by-step.



NOTICE

When the overheating of the drive unit occurs several times in a row, we recommend extending the duration of the cooling-phase significantly to avoid permanent damage to the device.

Whether and in what time the engine may overheat depends on the rated input and the present self-cooling power of the drive unit. The rated input of the system is determined by the following factors:

- Set speed
- Tool size and tool design
- Condition of the medium to be processed (such as: viscosity, temperature, volume, inclusions of air or solids, etc.)

Intensified conditions (e.g. higher ambient temperature, fluid temperature or differing viscosity levels) can shift the time of overheating of the drive unit. Therefore, careful handling of the X 1000 drive unit is fundamental!

9 Maintenance and Cleaning

9.1 Cleaning the Drive Unit UNIDRIVE X 1000/X 1000D

To clean the drive unit use only water with a detergent that contains tensides or use isopropyl alcohol for stubborn soiling.

9.2 Cleaning the Homogenising Tool



⚠ WARNING

The generators have sharp edges. Handle with care. Risk of injuries!

To avoid clogging of the shaft, clean shaft, generator and in case of G-shaft the sealing parts after each use. This is done by operating it in a solvent which dissolves substance residues and is not harmful to the gasket. This is usually sufficient to clean the Generator.

Chemical sterilization may be also a method. General-purpose disinfectants such as formalin, alcohol, etc. may be used. It is important to remove disinfectant residues with sterilized water.

NOTICE

Make sure that the bearings, O-rings and gaskets are resistant to solvents.

Sterilization with moist heat:

This denotes the use of a steam jet pressurized to 2 bars at 120°C.

9.3 Rough Cleaning of Flow-Through-Chamber and Tubing

Flow-through-chamber and tubing might be roughly cleaned by pumping a suitable liquid through the chamber instead of the medium.

9.4 Cleaning the Inside of the Flow-Through-Chamber

Before cleaning please disconnect first drive unit from mains and detach the flow through chamber DK 30

The flow-through chamber can be disassembled without disconnecting the tubes. **Attention:** Rest liquid escapes!!!

Please take care that neither O-rings nor any other parts experience any damage during assembling respectively disassembling the flow-through chamber.

Before operation check the slip ring sealing. The slip ring should be easily turned by hand. If it is stuck after a long period without use just loosen slip ring by wetting it with water. Never operate the flow-through chamber if the slip ring is stuck!!

NOTICE

Always unscrew the rotor before the stator.

Note: Rotor and stator wrenches are listed in the price list under accessories.

9.5 Maintenance of drive units UNIDRIVE X 1000/X 1000D

| | | |
|---|-----------------|--|
|  | ⚠WARNING | Do not open the instrument. Repairs are only to be carried out by trained service technicians. |
|---|-----------------|--|

9.6 Maintenance homogenizing tools

| | | |
|---|-----------------|--|
|  | ⚠WARNING | The generators have sharp edges. Handle with care. Risk of injuries! |
|---|-----------------|--|

The gaskets in the homogenizing tools must be constantly monitored. In the event of leakage the suction effect of the rotating Shaft can cause the medium to penetrate as far as the drive unit. If liquid emerges from the side hole at the top of the Shaft tube stop work immediately and check the gaskets. The function of the Generators depends on the condition of the sharp edges on the rotor and stator. These edges may be blunted very quickly in abrasive media, reducing the effectiveness of homogenizing.

10 Dismantling, Transport and Storage

10.1 Dismantling

1. Switch the unit off.
2. Disconnect it form the mains supply.
3. Remove any glass beakers and any other equipment around the drive unit.
4. Remove the homogenizing tool by loosening the holding screw (6).
5. Loosen the special clamp and remove the clamp from the support rod.
6. Unscrew the support rod.
7. Now you may remove the instrument from the working area.

10.2 Transport and Storage

Prior to transport:

Switch the instrument off and proceed with dismantling as described under "Dismantling".

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)
- Storage ambient:

Max. ambient temperature : RT to +40°C

Max. humidity: 80%

11 Disposal



Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: Metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

12 Warranty and Liability

The manufacturer agrees to either repair, or replace, at the manufacturer's discretion, any defects in materials or workmanship which develop within 24 months of the delivery of this product to the original user. In the event of replacement, the replacement unit will be guaranteed for the remainder of the original twenty-four (24) months period or ninety (90) days, whichever is longer.

If this product should require service, contact your local distributor or manufacturer for necessary instructions.

This guarantee will not apply if the defect or malfunction was caused by accident, neglect, unreasonable use or fitness for a particular purpose, which extend beyond the description and period set forth herein.

The manufacturer's sole obligation under this guarantee is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind, resulting from use or possession of the product.



WARNING

The user has to determine, if the instrument is suitable for his specific application. If there are any further queries, contact your local dealer or the manufacturer direct.

13 Technical Data

13.1 Technical Data UNIDRIVE X 1000/X 1000D

| Type: UNIDRIVE X 1000 | Specifications |
|---------------------------------|--------------------------|
| Motor Power intake | 1050 Watt |
| Motor Power output | 650 Watt |
| Idle speed | 4000 - 33000 rpm |
| Voltage | 230V/50 Hz or 115V/60 Hz |
| Case dimensions | 75x75x290 mm, WxDxH |
| Speed Display | - |
| Weight | 2.1 kg |
| Permissible ambient temperature | 5 - 40 °C |
| Permissible humidity | 80 % RH |

| Type: UNIDRIVE X 1000D | Specifications |
|---------------------------------|--------------------------|
| Motor Power intake | 1050 Watt |
| Motor Power output | 650 Watt |
| Idle speed | 4000 - 33000 rpm |
| Voltage | 230V/50 Hz or 115V/60 Hz |
| Case dimensions | 75x75 x 300 mm, WxDxH |
| Speed Display | Yes |
| Weight | 2.3 kg |
| Permissible ambient temperature | 5 - 40 °C |
| Permissible humidity | 80 % RH |

14 Repairs

| | |
|---|--|
|  | WARNING When returning instruments for repair that have come into contact with hazardous substances, please: Fill in attached "Repair Return Form" Provide precise information on the relevant medium Take protective measures to ensure the safety of our receiving and maintenance personnel Mark the package as appropriate for hazardous materials. |
|---|--|

14.1 Repair Return Form

| CONTACT/USER INFORMATION | | |
|---|--------------------|----------------|
| Contact: | Phone No. | |
| Fax No. | Email: | |
| Billing: | Shipping: | |
| Company | Company | |
| Address | Address | |
| INSTRUMENT INFORMATION | | |
| Model | Serial No. | |
| Please describe all problems/malfunctions | | |
| Operating Conditions (please fill in if applicable) | | |
| Ambient Temp. | Humidity | Speed |
| Load | Volume | Viscosity |
| Temperature in °C | Sample Temperature | Operating Time |
| Sample Description* | | |

*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to Ingenieurbüro CAT, M. Zipperer GmbH and an MSDS for hazardous material must be included with the instrument.

| RETURN SHIPPING | | | |
|------------------------------|--|------------------------------------|----------------------------------|
| UPS <input type="checkbox"/> | Air Parcel Post <input type="checkbox"/> | Collect** <input type="checkbox"/> | Other** <input type="checkbox"/> |

**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

| PACKAGING INSTRUCTIONS TO RETURN AN INSTRUMENT FOR REPAIR | |
|---|--|
|---|--|

- ✓ Remove all accessories (e.g. homogeniser tools, stirring paddles) from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenisers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.
Example: "Handle with care" or "Fragile- Delicate Instrument"

Ingenieurbüro CAT
M. Zipperer GmbH
Wettelbrunner Str. 6
D-79282 Ballrechten-Dottingen
Tel.: ++49-(0)7634-5056-800
Fax: ++49-(0)7634-5056-801
www.cat-ing.de
info@cat-ing.de

Gebrauchsanleitung

Dispergiersystem UNIDRIVE X 1000 und X 1000D

1 Allgemeine Information

Das Dispergiergerät Unidrive X 1000/Unidrive X 1000D wurde gem. DIN EN 61010 gefertigt und geprüft. Das Gerät hat das Werk in sicherheitstechnisch einwandfreiem Zustand verlassen. Um diesen Zustand zu erhalten und einen gefahrlosen Betrieb sicherzustellen, muss der Anwender die Hinweise und Warnvermerke beachten, die in dieser Gebrauchsanleitung enthalten sind.

2 Auspacken des Gerätes

Untersuchen Sie den Ummkarton. Achten Sie auf Beschädigungen wie z.B. eingedrückte Wände chemische Beeinträchtigungen des Kartons, Wassermarken, oder andere physikalische Einflüsse, die den Inhalt beschädigt haben könnten. Bei Beschädigungen informieren Sie bitte sofort den Spediteur und fragen Sie nach den notwendigen Maßnahmen. Wenn es keine sichtbaren Beschädigungen des Ummkartons gibt, dann öffnen Sie diesen vorsichtig.

2.1 Lieferumfang und Zubehör

Eine vollständige Sendung hat folgenden Umfang:

| Typ | Beschreibung | Art.Nr. |
|----------------|-------------------------|-----------------|
| X 1000 | Dispergierantrieb, 230V | 60402-0050 oder |
| X 1000 | Dispergierantrieb, 115V | 60402-0051 oder |
| X 1000D | Dispergierantrieb, 230V | 60406-0050 oder |
| X 1000D | Dispergierantrieb, 115V | 60406-0051 oder |
| | 1 Haltestange | |
| | 1 Gebrauchsanleitung | |

2.2 Dispergierwerkzeuge und Zubehör (separat zu bestellen):

| T 6 | T 6 long | T 10 | T 17 | T 20 | G 20 |
|---|---|--|---|---|------------|
|  |  |  |  |  | |
| 2,0 - 50 ml | 10,0 - 250 ml | 250 - 1500 ml | 250 - 1500 ml | 250 - 5000 ml | |
| 55 mm | 120 mm | 150 mm | 150 mm | 200 mm | 215 mm |
| fine F | | | | | |
|  |  |  |  |  | |
| 60420-000F | 60410-000F | 60421-000F | 60426-000F | 60422-000F | 60423-000F |
| low viscosity N | | | | | |
|  |  |  |  |  | |
| 60420-000N | 60410-000N | 60421-000N | 60426-000N | 60422-000N | 60423-000N |
| viscous V | | | | | |
|  |  |  |  |  | |
| 60420-000V | 60410-000V | 60421-000V | 60426-000V | 60422-000V | 60423-000V |
| knife M | | | | | |
|  |  |  |  |  | |
| 60420-000M | 60410-000M | 60421-000M | 60426-000M | 60422-000M | 60423-000M |

T = Schaft mit PTFE Lager: Dieser Schaft-Typ ist für Standard Anwendungen. (für flüssige Medien wie z.B. Wasser, Öl, o.ä.) Das Lager muss von dem Medium, welches in das Lager eindringt, gekühlt werden. Ein Trockenlauf ist nicht zulässig.

G = Schaft mit Gleitringdichtung: Dieser Schaft Typ ist sehr gut für aggressive oder abrasive Medien geeignet oder wenn Druck aufgebaut wird (z.B. bei viskosen Medien) oder wenn Feststoffe zerkleinert werden sollen.
Das Lager ist abgedichtet und für einen Trockenlauf geeignet

| Werkzeuge | Beschreibung | Art Nr. |
|--------------------------|--|------------|
| Plattenstativ inkl. Stab | Plattformabmessung: 330x200 mm, Stab-Durchmesser: 16mm, Länge:690mm | 60491-0000 |
| Kreuzklemme | Halteklemmer zur Befestigung am Stativ Material: Aluminiumguss | 60492-0000 |
| Kreuzklemme | Halteklemmer zur Befestigung am Stativ Material: Stahl | 60492-0050 |
| Universal-Schlüssel | Universal-Rotor-Schlüssel für 6 - 30 mm | 60470-0000 |
| Steckschlüssel | für alle Schäfte | 60471-0000 |
| Rotorschlüssel | für Schaft T 10 | 60472-0000 |



Wichtiger Hinweis:

Bitte lesen Sie die Gebrauchsanleitung gründlich durch und stellen Sie sicher, dass Sie alles verstanden haben, bevor Sie das Gerät in Betrieb nehmen. Sollten Sie nach dem Durchlesen der Gebrauchsanleitung noch weitere Fragen zu Installation, Betrieb oder Wartung haben, wenden Sie sich bitte an den Lieferanten oder den Hersteller des Gerätes unter folgender Adresse:

Ingenieurbüro CAT
M. Zipperer GmbH
Wettelrunner Str. 6
D-79282 Ballrechten-Dottingen
Tel.: ++49-(0)7634-5056-800
Fax: ++49-(0)7634-5056-801
www.cat-ing.de
info@cat-ing.de



Wichtiger Hinweis:

Diese Gebrauchsanleitung sollte jederzeit zur Verfügung stehen, insbesondere demjenigen, der dieses Gerät nutzen möchte. Deshalb sollte diese Gebrauchsanleitung in der Nähe des Gerätes aufbewahrt werden.



Wichtiger Hinweis:

Stellen Sie sicher, dass das Gerät nur gemäß den zulässigen Bedingungen beschrieben unter dem Punkt "Technische Daten" betrieben wird.



Die Geräte sind standardmäßig mit einem Eurostecker EN 50075:1990, 2,5 A, 250 V versehen. Für Nordamerika mit einem US-Norm Stecker (NEMA 1 ,Typ A) 10 A, 125 V.



Falls Sie das Gerät in einem Land mit anderem Stecker-System betreiben möchten, müssen Sie einen zugelassenen Adapter verwenden oder der mitgelieferte Stecker muss durch einen Fachmann ausgetauscht und durch einen für dieses Netz passenden und zugelassenen Stecker ersetzt werden.

3 Sicherheitshinweise

3.1 Erläuterung der Sicherheitshinweise in der Gebrauchsanleitung



Dieses Symbol hebt folgende Hinweise hervor, welche durch den Anwender unbedingt beachtet werden müssen! Jegliche Missachtung der nebenstehenden Hinweise, kann die sichere Funktion wie die Sicherheit des Anwenders selbst beeinträchtigen.



Dieses Symbol hebt Verbote hervor, welche durch den Anwender unbedingt beachtet werden müssen! Jegliche Missachtung der nebenstehenden Verbote, kann die Funktionstüchtigkeit, sowie die Sicherheit des Anwenders erheblich beeinträchtigen.



Dieses Symbol hebt Hinweise hervor, welche durch den Anwender unbedingt beachtet werden sollten um einen sicheren Betrieb des Gerätes gewährleisten zu können.



Achtung:
Brand- oder Explosionsgefahr!

3.2 Erläuterung der Sicherheitshinweise auf dem Gerät

Warnung vor allg. Gefahren:



Dieses Symbol weist darauf hin, dass die Gebrauchsanleitung unbedingt vor Inbetriebnahme des Gerätes vollständig gelesen und verstanden werden sollte. Eine unsachgemäße Handhabung kann den sicheren Betrieb des Gerätes und die Gesundheit des Anwenders gefährden.



Dieses Symbol weist darauf hin, dass es unbedingt notwendig ist Augen- und Hörschutz zu tragen, um Gesundheitsschäden zu vermeiden.

3.3 Allgemeine Sicherheitshinweise



Achtung:

Beachten Sie alle im Labor geltenden Sicherheits- und Unfallverhütungsvorschriften!



Achtung:

Verwenden Sie dieses Gerät nicht in unsicherer Umgebung und speziell nicht in explosionsgefährdender Umgebung. Lebensgefahr!



Achtung:

Ausschließlich unterwiesene Anwender dürfen das Gerät in Betrieb nehmen.



Achtung:

Achten Sie vor der Verbindung des Gerätes mit dem Stromnetz darauf, dass die Netzspannung mit den Angaben auf dem Typenschild übereinstimmt.



Achtung:

Äußerste Vorsicht beim Umgang mit leicht entzündlichen Medien. Beachten Sie die Sicherheitsdatenblätter. Das Gerät ist nicht explosionsgeschützt!



Der Hauptschalter trennt das Gerät nicht vollständig von der Stromquelle. Um das Gerät vollständig vom Netz/der Stromversorgung zu trennen, ziehen Sie bitte den Netzstecker.



Hinweis:

Stellen Sie sicher, dass alle Sicherheitshinweise des Gerätes bzw. auf dem Gerät selbst, während des Betriebs deutlich zu sehen sind.



Achtung:

Öffnen Sie das Gerät nicht. Reparaturen sind eingewiesenen Service-Technikern vorbehalten.



Bitte achten Sie auf sichere Befestigung des Gerätes.



Tragen Sie einen Gehörschutz beim Überschreiten des Geräuschpegels von 85dB(A).



Schützen Sie während des Arbeitens Ihrer Augen mit einer geeigneten Schutzbrille.



Stellen Sie sicher, dass das Gerät an einem Stativ befestigt ist!

4 Aufbau und Betrieb


4.1 Zweckbestimmung

Der Dispergierantrieb UNIDRIVE X 1000/X 1000D und geeignete von uns empfohlene Dispergierwerkzeuge bilden ein System zum Bearbeiten von fließfähigen bzw. flüssigen Medien im Chargenbetrieb. Das System wird zum Dispergieren sowie Emulgieren verwendet.

Verwendungsgebiete sind Laboratorien sowie Bildungseinrichtungen.


Für einen sicheren Betrieb muss das System an einem geeigneten Stativ befestigt werden.

4.2 Beschreibung des Antriebs UNIDRIVE X 1000/X 1000D



4.3 Aufbau

4.3.1 Montage der Haltestange




- Schrauben Sie die Haltestange (8) in die Gewindebohrung (7) auf der Rückseite des Gerätes.
- Stecken Sie einen geeigneten Schraubendreher in die Bohrung am Ende der Haltestange drehen Sie nach rechts und ziehen somit die Haltestange fest.

4.3.2 Befestigen des Antriebs am Stativ




Um einen sicheren Betrieb zu gewährleisten, muss der Antrieb mit Hilfe einer Kreuzklemme (9) am Stativ befestigt werden. Vibrationen können dazu führen, dass sich die Schrauben (10) an der Klammer lösen. Überprüfen Sie regelmäßig vor jedem Gebrauch, ob die Schrauben der Kreuzklemme noch fest angezogen sind.



- Befestigen Sie die Kreuzklemme (9) am Stativstab.
- Jetzt befestigen Sie den Antrieb am Stativ indem Sie die Haltestange (8) in die freie Öffnung der Kreuzklemme (9) schieben.
- Drehen Sie die Schrauben (10) der Kreuzklemme (9) fest zu.

4.4 Beschreibung eines Dispergierwerkzeugs




4.4.1 Anbringen des Dispergierwerkzeugs



Die Generatoren haben scharfe Kanten. Verletzungsgefahr!


Dispergierwerkzeuge werden immer zusammengebaut angeliefert. Der Antriebsstift (11) befindet sich am oberen Ende des Schaftes. Er dient zur Verbindung des Schaftes mit dem Antrieb. Führen Sie den Schaft in die Kupplung des Antriebs ein. Jetzt drehen Sie den Schaft etwas hin und her bis Sie den Antriebsstift im Kupplungsstück einrasten hören. Befestigen Sie jetzt den Schaft mit der Befestigungsschraube (6).



Lassen Sie das Dispergierwerkzeug nie trocken laufen! Das zu bearbeitende Medium dient gleichzeitig auch zur Kühlung der Lager und Dichtungen. Hierdurch wird eine Beschädigung des Werkzeugs vermieden. Der Abstand zwischen Dispergierwerkzeug und Gefäßboden darf nicht weniger als 10mm betragen. Der Flüssigkeitsstand sollte mindestens 55 mm betragen. Um einen Vortex zu vermeiden muss das Dispergierwerkzeug außer-mittig ins Gefäß eingetaucht werden. Sind alle diese Voraussetzungen erfüllt ist das System betriebsbereit.

4.5 Betrieb UNIDRIVE X 1000/X 1000D

| | |
|---|--|
|  | Achtung: Achten Sie vor der Verbindung des Gerätes mit dem Stromnetz darauf, dass die Netzspannung mit den Angaben auf dem Typenschild übereinstimmt. |
|  | Achtung: Das Gerät hat einen internen Thermischen Überlastschutz. Nach dem Ansprechen dieses Schutzes mindestens 10 Min. bis zum Wiedereinschalten warten, damit der Motor abkühlen kann. |



4.5.1 Anschließen des Gerätes ans Stromnetz

Wenn das System komplett zusammengebaut und am Stativ befestigt ist, stecken Sie den Stecker in eine geeignete Steckdose.

Das Gerät wird mit dem Ein-/Ausschalter (3) vorne am Gehäuse gestartet.

Position I: Motor läuft

Position O: Motor aus

Stellen Sie die Drehzahl am Drehzahlsteller (2) ein. Starten Sie immer mit der Mindestdrehzahl und erhöhen Sie diese langsam bis zur gewünschten Einstellung.



Achtung:

Der Ein-/Ausschalter (3) ist nur einpolig und trennt das Gerät nicht vollständig vom Netz. Zur vollständigen Trennung vom Netz, Stecker ziehen!



Bevor Sie das Gerät einschalten stellen Sie sicher dass:

- die Mindestdrehzahl eingestellt ist
- der Generator (12) des Dispergierwerkzeugs in der Flüssigkeit eingetaucht ist. Um einen Vortex zu vermeiden, tauchen Sie den Schaft außer-mittig ein
- das Gefäß gegen Eigendrehung, Schütteln und Wandern gesichert ist.
- Die Kühlslitze nicht abgedeckt sind, da hiermit der Antrieb gekühlt wird.

4.5.2 Beschreibung der Stellradzahlen des X 1000 und X 1000D:

| | | |
|---|---|-------------|
| 1 | = | 4000 1/min |
| 2 | = | 8500 1/min |
| 3 | = | 15000 1/min |
| 4 | = | 20000 1/min |
| 5 | = | 25000 1/min |
| 6 | = | 33000 1/min |

Die eingestellte Drehzahl des UNIDRIVE X 1000D kann auf der Drehzahlanzeige (5) abgelesen werden.

4.6 Hinweise zum Dauerbetrieb der X 1000-Dispergierantriebe

Der X 1000-Dispergierantrieb ist bei Verwendung der Schäfte T6, T10, T17 und T/G20 für den Dauerbetrieb ausgelegt.

4.7 Überhitzungsschutz der X 1000-Dispergierantriebe

Zum Schutz vor Überhitzung verfügt der X 1000-Dispergierantrieb eine Sicherheitselektronik, welche bei Erreichen einer kritischen Temperatur in den Abkühlmodus (temperaturabhängiger Überlastschutz) schaltet. Die Maschine läuft dann mit stark reduzierter Drehzahl weiter und die Konstant-Elektronik wird deaktiviert. Nach einer Abkühlzeit von 60 Sekunden ist die Maschine wieder betriebsbereit. Hierfür die Maschine aus und einschalten, um die Konstant-Elektronik zu aktivieren. Vor dem Einschalten über das Stellrad die Mindestdrehzahl einstellen um ein Verspritzen des Mediums zu vermeiden. Die Drehzahl kann dann per Stellrad wieder auf die gewünschte Geschwindigkeit kontinuierlich hochgefahren werden.



Wenn eine Überhitzung des Motors von Lauf zu Lauf immer frühzeitiger eintritt, empfehlen wir die Dauer der Abkühlzeit deutlich zu erhöhen um dauerhaften Schaden am Gerät zu vermeiden.

Ob und in welcher Zeit sich der Motor überhitzen kann hängt von der dauerhaften Leistungsaufnahme und der vorliegenden Eigenkühlleistung des Motors ab. Die Leistungsaufnahme des Systems wird durch folgende Faktoren bestimmt:

- eingestellte Drehzahl
- Größe und Ausführung des Werkzeugs
- Beschaffenheit des zu bearbeitenden Mediums (wie z.B.: Viskosität, Temperatur, Volumen, Einschlüsse von Luft oder Festkörpern usw.)

D.h.: Liegen erschwerte Bedingungen (z.B.: höhere Umgebungs-, Medientemperatur oder abweichende Medienviskosität) vor, kann sich eine Überhitzung des Motors früher oder später einstellen. Daher ist eine umsichtige Handhabung der X 1000-Dispergierantriebe eine grundlegende Voraussetzung!

4.8 Austausch der Dispergierwerkzeuge



Achtung die Generatoren haben scharfe Kanten. Verletzungsgefahr!

- Bevor Sie das Dispergierwerkzeug wechseln, schalten Sie das Gerät aus und trennen es vom Netz, indem Sie den Stecker ziehen.
- Befestigungsschraube (6) lösen und Schaft mit der Hand aus der Kupplung ziehen.
- Führen Sie den neuen Schaft in die Kupplung des Antriebs ein. Jetzt drehen Sie den Schaft etwas hin und her bis Sie den Antriebstift im Kupplungsstück einrasten hören.
- Befestigen Sie jetzt den Schaft mit der Befestigungsschraube (6).

4.8.1 Austausch des Generators



Die Generatoren haben scharfe Kanten. Achtung: Verletzungsgefahr!

Um den Generator auszutauschen (Rotor und Stator), entfernen Sie den Schaft aus dem Antrieb wie unter Punkt 4.6 beschrieben.

Den Rotor lösen Sie mit Hilfe der zwei Spezialschlüssel (separat zu bestellen). Der Steckschlüssel dient zum Gegenhalten des Schaftes, der Rotorschlüssel zum Lösen bzw. Anziehen.

- Stecken Sie den Rotorschlüssel in den Rotor und den Antriebstift in den Steckschlüssel.
- Drehen Sie den Rotorschlüssel nach links, um den Rotor zu lösen.
- Danach kann der Stator von Hand durch Drehen nach **rechts** gelöst werden (Achtung: Linksgewinde!)



Wichtiger Hinweis

Lösen Sie den Rotor immer zuerst bevor der Stator abgeschraubt wird.

5 Reinigung und Wartung

5.1 Reinigung des Antriebs UNIDRIVE X 1000/X 1000D

Zur Reinigung des Antriebs verwenden Sie bitte nur Wasser und ein geeignetes Reinigungsmittel (z.B. Geschirrspülmittel) oder benutzen Sie Isopropylalkohol um hartnäckige Verschmutzungen zu entfernen.

5.2 Reinigen des Dispergierwerkzeugs



Die Generatoren haben scharfe Kanten. Vorsicht! Verletzungsgefahr.

Um ein Verstopfen zu vermeiden, reinigen Sie Schaft, Generator und bei Verwendung eines G-Schafts die Dichtungsteile nach jedem Gebrauch. Hierzu tauchen Sie den Schaft in ein geeignetes Lösungsmittel, welches das Werkzeug von Rückständen befreit und das die Dichtungen nicht angreift.

Dies ist normalerweise ausreichend um den Generator zu reinigen.

Chemische Sterilisation kann zur Reinigung ebenfalls angewendet werden. Allgemein gebräuchliche Desinfektionsmittel wie z.B. Alkohol können hierzu verwendet werden. Wichtig: Entfernen Sie Desinfektionsmittel-Rückstände mit sterilisiertem Wasser.



Wichtiger Hinweis:

Stellen Sie sicher, dass Lager, O-Ringe und Dichtungen chemisch-beständig zu den verwendeten Lösungsmitteln sind.

Dampfsterilisation:

Benutzen Sie hierzu einen Dampfstrahl mit 2 bar Druck bei 120°C.

5.3 Vorreinigung der Durchflusskammer und den Schläuchen

Durchflusskammer und Schläuche können durch pumpen einer geeigneten Flüssigkeit vorgereinigt werden.

5.4 Reinigung des Kammerinnenraums

Vor der Reinigung entfernen Sie die Durchlaufkammer DK 30 vom Dispergiersystem.

Die Durchlaufkammer kann vom Dispergiersystem entfernt werden ohne dass vorher die Schläuche abgenommen werden. **Achtung:** Restflüssigkeit kann ausfließen.

Bitte achten Sie darauf, dass keine Teile während der Montage beschädigt werden.

Achtung: Überprüfen Sie vor Gebrauch die Gleitringdichtung. Die Dichtung sollte sich leicht von Hand drehen lassen. Sitzt diese fest, lösen Sie diese, indem Sie sie mit Wasser befeuchten. Betreiben Sie die Durchflusskammer nicht wenn die Gleitringdichtung des Schaftes festsitzt.

Hinweis: Rotor- und Steckschlüssel sind als Zubehör in der Preisliste aufgeführt.

5.5 Wartung der Dispergierwerkzeuge



Die Generatoren haben scharfe Kanten. Achtung! Verletzungsgefahr!

Die Dichtungen der Dispergierwerkzeuge müssen ständig überwacht werden. Bei Leckage kann der Ansaugeffekt des rotierenden Schaftes bewirken, dass Flüssigkeit bis in den Antrieb gelangt.

Wenn Flüssigkeit aus den Seitenlöchern am oberen Ende des Schafts heraustritt, stoppen Sie den Betrieb sofort und überprüfen Sie die Dichtungen. Die Funktion des Generators hängt von dem Zustand der Zähne des Generators ab. Diese Zähne können beim Bearbeiten von abrasiven Medien sehr schnell verschleiß. Dies vermindert die Effektivität des Vorgangs.

6 Abbau, Transport und Lagerung

6.1 Abbau

1. Schalten Sie das Gerät aus.
2. Ziehen Sie den Netzstecker.
3. Entfernen Sie sämtliche Behälter und Gegenstände um das Gerät herum.
4. Entfernen Sie das Dispergierwerkzeug aus dem Antrieb.
5. Lösen Sie die Schraube der Kreuzklemme und nehmen den Antrieb vom Stativ.
6. Entfernen Sie die Kreuzklemme von der Haltestange.
7. Schrauben Sie die Haltestange heraus.
8. Jetzt können Sie das Gerät aus der Arbeitsumgebung entfernen.

6.2 Transport und Lagerung

Vor dem Transport:

- Schalten Sie das Gerät aus und gehen Sie wie unter Punkt 8.1 beschrieben vor.
- Verpacken Sie Gerät und Zubehör in den Originalkarton oder einen anderen geeigneten Karton, um das Gerät während des Transports zu schützen.
- Verschließen Sie den Karton mit Klebeband.

Lagern Sie das Gerät in einer trockenen Umgebung. Beachten Sie bitte die angegebenen Umgebungsbedingungen wie Temperatur und Luftfeuchtigkeit.

Setzen Sie das Gerät keinen mechanischen Erschütterungen oder Vibrationen während des Transports aus.

Sollte der Originalkarton nicht mehr zur Verfügung sein markieren Sie die Verpackung wie folgt:

- Glassymbol (Vorsicht zerbrechlich)
- Regenschirm (trocken aufbewahren)
- Inhalt (Aufstellung des Inhalts)
- Umgebungsbedingung:
 - Max. Umgebungstemperatur : RT bis +40°C
 - Max. Feuchtigkeit: 80%

7 Entsorgung



Schadhafte und/oder entsorgte elektronische oder elektronische Geräte müssen an den dafür vorgesehenen Recycling-Stellen abgegeben werden.
Auch das Verpackungsmaterial sollte umweltgerecht (Materialtrennung) entsorgt werden.

8 Garantie und Haftungsausschluss

Der Hersteller verpflichtet sich dem Käufer gegenüber dieses Produkt entweder durch Reparatur oder nach Wahl des Herstellers durch Austausch jeden Fehler in Material oder Verarbeitung zu beheben, der sich innerhalb von 24 Monaten nach Auslieferung dieses Produktes herausstellt. Für den Fall des Austausches wird die Garantie weiterhin bis zum Ende der 24monatigen Laufzeit ab Kaufdatum laufen.

Für allfälligen Kundendienst ist der Händler oder der Hersteller einzuschalten.

Diese Garantie wird nicht gewährt, wenn der Defekt oder die Fehlfunktion verursacht wird durch Unfall, Nachlässigkeit, unsachgemäßen Gebrauch, unsachgemäßen Service und andere Gründe, die nicht auf fehlerhafte Materialien oder Verarbeitung zurückzuführen sind.

Sämtliche Garantie- und Gewährleistungsansprüche sind beschränkt auf die fachgerechte und sachgemäße Anwendung des Gerätes unter Beachtung der allgemeinen Sicherheitsvorschriften und dieser Gebrauchsanleitung. Der Hersteller ist zu keinerlei anderen Leistungen als den in dieser Garantieerklärung enthaltenen verpflichtet. Insbesondere sind Personen- und Sachschäden, die als Folge von fehlerhaften Teilen oder sonstigen Funktionsfehlern auftreten, von der Haftung ausgeschlossen.

Die Verpflichtung des Herstellers ist im Umfang begrenzt auf die Reparatur oder den Austausch von fehlerhaften Teilen. Der Hersteller ist unter keinen Umständen zu Schadensersatz verpflichtet für Folgeschäden jeder Art, die aus der Verwendung und dem Besitz dieses Produktes sich ergeben können.

9 Reparaturen



Bitte sorgen Sie bei der Anlieferung von Reparaturgeräten, die mit gefährlichen Arbeitsstoffen in Berührung gekommen sind dafür:

- *Möglichst genaue Stoffangabe des Mediums*
- *Schutzmaßnahmen zum sicheren Umgang für unser Annahme- und Wartungspersonal*
- *Kennzeichnung der Verpackung gemäß der Gefahrenstoffverordnung*

10 Technische Daten

| Type: UNIDRIVE X 1000 | Spezifikationen |
|-------------------------------|----------------------------|
| Motor Leistungsaufnahme | 1050 Watt |
| Motor Leistungsabgabe | 650 Watt |
| Leerlauf | 4000 - 33000 1/min |
| Netzspannung | 230V/50 Hz oder 115V/60 Hz |
| Abmessungen | 75x75x290 mm, BxTxH |
| Drehzahlanzeige | - |
| Gewicht | 2.1 kg |
| Zulässige Umgebungstemperatur | 5 - 40 °C |
| Zulässige Luftfeuchtigkeit | 80 % RH |

| Type: UNIDRIVE X 1000D | Spezifikationen |
|-------------------------------|-----------------------------|
| Motor Leistungsaufnahme | 1050 Watt |
| Motor Leistungsabgabe | 650 Watt |
| Leerlauf | 4000 - 33000 Upm |
| Netzspannung | 230V/50 Hz oder 115V/60 Hz |
| Abmessungen | 75 x 75 x 300 mm, B x T x H |
| Drehzahlanzeige | Ja |
| Gewicht | 2,3 kg |
| Zulässige Umgebungstemperatur | 5 - 40 °C |
| Zulässige Luftfeuchtigkeit | 80 % RH |

10.1 Rücksendeformular

Wichtig:

1. Kalibrieren des Instruments wird automatisch durchgeführt.
2. Sie erhalten einen Kostenvoranschlag für die zu erwartenden Reparaturkosten.
3. Reparaturen werden bei Freigabe durchgeführt.

KONTAKT

| | |
|---------------------|------------------|
| Ansprechpartner: | Tel. |
| Fax Nr. | Email: |
| | |
| Rechnungsanschrift: | Lieferanschrift: |
| Firma | Firma |
| Adresse | Adresse |
| | |
| | |
| | |
| | |

GERÄTEINFORMATION

| | | |
|-----------------------------|-------------|--|
| Modell | Serien- Nr. | |
| Problem-/Fehlerbeschreibung | | |
| | | |
| | | |

Betriebsbedingungen (bitte ausfüllen, wenn zutreffend)

| | | |
|---------------------|-------------------|------------|
| Umgebungstemp. | Rel. Feuchtigkeit | Drehzahl |
| Zuladegewicht | Volumen | Viskosität |
| Temperatur in °C | Probentemperatur | Laufzeit |
| Probenbeschreibung* | | |

*HINWEIS: Sollte das Gerät mit gefährlichen Stoffen in Berührung gekommen sein, muss es vor dem Versand an Ingenieurbüro CAT, M. Zipperer GmbH dekontaminiert werden. Fügen Sie der Rücklieferung die entsprechenden Sicherheitsdatenblätter bei.

RÜCKSENDUNG

| | | |
|------------------------------|------------------------------------|-----------------------------------|
| UPS <input type="checkbox"/> | Paketpost <input type="checkbox"/> | Andere** <input type="checkbox"/> |
|------------------------------|------------------------------------|-----------------------------------|

**Bitte nennen Sie uns Ihren Spediteur, Anschrift, Sachbearbeiter/in sowie dessen/deren Telefonnummer

VERPACKUNGSHINWEISE

- ✓ Bitte entfernen Sie sämtliches Zubehör. (Z.B. Dispergierwerkzeuge, Rührpaddel, etc.)
- ✓ Reinigen Sie das Gerät und Zubehör von Probematerial-Rückständen
- ✓ Legen Sie dem Gerät die entsprechenden Sicherheitsdatenblätter für eventuell benutzte gefährliche Stoffe bei.
- ✓ Verpacken Sie das Gerät in die Originalverpackung. Sollte diese nicht mehr vorhanden sein, achten Sie darauf, dass das Gerät sowie Zubehör ausreichend geschützt verpackt wird.
- ✓ Plattenstativ und Halteklammer bitte nicht einschicken. Diese senden Sie nur zurück, wenn sie beschädigt sind. Um Transportschäden zu vermeiden, achten Sie bitte darauf, dass Gerät und Zubehör individuell verpackt werden.
- ✓ Verpacken Sie das Gerät mit dem Zubehör in einen ausreichend starken Versandkarton. Kennzeichnen Sie dessen Außenseite mit
Beispiel: "Vorsicht" oder "Vorsicht-zerbrechlich" und schicken Sie die Waren an:

Ingenieurbüro CAT
M. Zipperer GmbH
Service Department
Wettelbrunner Str. 6
D-79282 Ballrechten-Dottingen



EG - Konformitätserklärung nach Maschinenrichtlinie 2006/42/EG

EC - Declaration of conformity in accordance with the EEC machine directive 2006/42/EG appendix II 1A

Déclaration de conformité CE selon Directive Européenne 2006/42/EG relative aux machines appendix

II 1A

Wir erklären in alleiniger Verantwortung, dass dieses Produkt (siehe Tabelle) den Bestimmungen der Richtlinien entspricht und mit den folgenden Normen oder normativen Dokumenten übereinstimmt:

We declare under our sole responsibility that this product (see table) corresponds to the regulations and conforms with the standards or standardized documents:

Nous déclarons sous notre propre responsabilité que ce produit est en conformité avec les normes ou documents normalisés suivant:

Folgende weitere EU-Richtlinien wurden angewandt :

The agreement with further valid guidelines/regulations following for the product is explained:

Les autres directives européennes suivantes ont été appliquées:

EMV – Richtlinie 2014/30/EG

RoHS – Richtlinie 2011/65/EU

WEEE II – Richtlinie 2012/19/EU

Folgende harmonisierte Normen wurden angewandt :

Applied one harmonized standards in particular:

Les normes harmonisées suivantes ont été appliquées:

DIN EN 61326-1 :2013-07 DIN EN ISO 12100 :2011-03

DIN EN 61010-1:2011-07 DIN EN 61010-2-051 :2016-02

Bei einer nicht mit uns abgestimmten Änderung des Gerätes verliert diese Erklärung ihre Gültigkeit und die Gewährleistung erlischt.

In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

Dans le cas d'une modification de l'appareil qui n'a pas été convenue avec nous, cette déclaration perd sa validité et la garantie expire.

Ingenieurbüro CAT M. Zipperer GmbH


Alexander Schecklein
CE-Bevollmächtigter
79282 Ballrechten-Dottingen, den 08.08.2018

| Type | 230 Volt Version | | | 115 Volt Version | | |
|----------------|------------------|----------|--------|------------------|----------|--------|
| X 120 | 230 V + 10% | 50/60 Hz | 140 W | 115 V + 10% | 50/60 Hz | 140 W |
| X 360 | 230 V + 10% | 50/60 Hz | 350 W | 115 V + 10% | 50/60 Hz | 350 W |
| X 1000, X1000D | 230 V + 10% | 50/60 Hz | 1050 W | 115 V + 10% | 50/60 Hz | 1050 W |
| X 1740 | 230 V + 10% | 50/60 Hz | 1800 W | 115 V + 10% | 50/60 Hz | 1800 W |